

# Molift Air 200

The compact, powerful hoist optimized for long term care

oetac<sup>®</sup>



# Molift Air 200

Molift Air 200 is an affordable, well-designed and highly functional ceiling hoist for users up to 205 kg. Combined with Molift Rail System and Molift slings, this hoist solves most day-to-day hoisting situations in health care and other patient handling environments.

# **Comfort and quality**

Molift Air 200 combines the need for maximum comfort with efficient and safe hoisting situations. The design and user-friendliness bring an exceptional freedom, which enables a single caregiver to perform daily lifting situations, allowing more time for caring for the individual user. For the caregiver, a better working environment is possible without the need for manual lifting. The user experiences a dignified, safe and comfortable lifting.

# The flexible hoist

Molift Air 200 is suitable for sitting and horizontal transfers, as well as standing and gait training situations.

Molift Air 200 is used with fixed ceiling rails but can also be used for portable free-standing systems. It is easy to handle during installation and servicing, due to a very low weight and the smart quick release function.

Molift Air 200 offers flexible solutions such as remote control and possibilities for room-to-room transfer. Add to this the safety features, the large hoisting interval and the multiple range of sling bars. This hoist is the complete solution in all health care segments.



Product model	Includii	Including	
Molift Air 200 SLA Battery	Hoist, han	dcontrol, charger	26200
Molift Air 200 Litium Battery	Hoist, han	dcontrol, charger	26201
Facts			
Maximum user weight (SWL)			205 kg
Total weight			8.6 kg
Length			360 mm
Width			190 mm
Height			195 mm
Battery SLA		1	2 V SLA 2,9 Al
Battery Litium		25,6 V I	iFePO4 2,4 Ah
Battery charger		100-240 V	AC, 50-60 Hz
Battery charge time			6 hrs
Motor			24 V D0
Hoisting interval			3000 mm
Hoisting speed 75 kg			60 mm/sed

Facts	
Expected lifetime	30 000 lifts SWL/10 years
Molift Service Software	Included
Emergency lowering	Mechanical and electrical
Emergency stop	Reachable from floor with string. Reset possible from floor
Protection class; hand control & hoist motor	IP24





Clean / disinfect the product using a standard non-abrasive cleaning agent with a pH between 5 and 9, or with a 70% disinfectant solution.



# **Features**



#### Small hoist - great lifting capacity

Molift Air 200 is very small and weights only 8,6 kg. Nevertheless it has a lifting capacity (SWL) of 205 kg. The size and weight makes it very easy to operate and move for the carer during different transfer situations.



### Silent and soft operation

Molift Air 200 works very quietly and has a soft start and stop operation. This results in a comfortable and safe feeling for the user during the transfer situation.



#### Quick release device for trolley

Molift Air 200 is easy to handle during installation and servicing, due to a very low weight and the smart quick release function. With a simple grip, the motor is disconnected from its trolley and removed.



#### Intuitive hand control

Molift Air 200 has a new intuitive handcontrol with attachment possibility on sling bar and optional docking station for wall.

# **Accessories**



### 4-point Sling bar

Light weight aluminium, SWL 350 kg  $\,$ 

4-pt Sling bar 350 (S)	2530210
4-pt Sling bar 450 (M)	2530211
4-pt Sling bar 550 (L)	2530212
4-pt Sling bar 650 (XL)	2530213



#### **Molift Air Scale**

This scale can calculate BMI (Body Mass Index) with a simple touch. Approved Class III.

Scale SWL 300 kg......1840000



#### 2-point Sling bar

Light weight aluminium, SWL 350 kg

2-pt Sling bar 350	(S)	2530115
2-pt Sling bar 450	(M)	2530120
2-nt Sling har 550	(1)	2530125



#### Molift Stretcher MRI

Flexible scoop stretcher, ideal for horizontal hoisting situations at hospitals and institutions. SWL 300 kg

Molift Stretcher MRI.....2150108



# 2-point Pigtail Sling bar

Light weight aluminium, SWL 350 kg  $\,$ 

2-pt Pigtail Sling bar 350 (S)	1830033
2-pt Pigtail Sling bar 450 (M)	1830032
2-pt Pigtail Sling bar 550 (L)	1830031



# Other accessories

Room-to-room unit ......2510295

3



# 8-point Sling bar

For transfer situations with a fabric stretcher. Steel, SWL 300 kg

8-pt Sling bar.....2140003

www.etac.com Molift Air 200



Etac is a world-leading developer and provider of ergonomic assistive devices and patient handling equipment. Our heart lies in the solutions that optimise quality of life for the individual, their family and caregivers.

For the latest news and continuously updated product information – visit: www.etac.com



