

Rose Comfort Max tilt and recline chair with casters

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D




FIGURE E




FIGURE F



FIGURE G



FIGURE H

- 38 degrees of tilt- lever on the back left positioned easily for the caregiver to adjust
- Recline up to 90 degrees (trendelenburg). Recline lever is Red and located on the back Right
- Seat Width 20" and seat depth 20"
- Seat height adjusts from 17"-22" in 1" increments
- Leg rest: Elevating, length adjustable, lever in the back or side positioned easily for caregiver to adjust
- Footrest: Removable single flip up footrest with angle adjustable foot plate, need also secure so it does not get lost
- Wings: Adjustable Removable for Upper Lateral Support
- Casters: (2) Front 5 total lock (12.5 cm) Rear (1) 5 Directional locking & (1) 5 total lock caster (12.5 cm)
- Standard Padding: Headrest cushion, shoulder bolsters, lower side pads, upper side pads, seat and back pad with cover, calf pad, sole pad
-
-

PRODUCTS

Description

CLN400-20	Chair with Caster Full Tilt and Recliner
-----------	--

SPECIFICATIONS

Description

Back Height	33"
Color	Blue
Dimensions	Arm Pad: 14" (L) x 3" (W) Leg Pad: 19" (L) x 12.5" (H) Foot Pad: 12" (W) x 22" (H)
Overall Height	23.5"
Overall Width	65"
Recline	Up to 90 degrees (trendelenburg)
Seat Depth	18"
Seat Height	Standard position: 17" - 22" In full recline: 21.5"
Seat Width	20"
Tilt Angle Range	Up to 38 degrees
Warranty	Lifetime (10 years on frame)
Weight Capacity	400 lbs
Width	Between Shoulder Width: 22.5"
Material Composition	16 guage tubular steel frame, 1 1/2" vinyl straps, flat head screws, foam

ACCESSORIES

Description

CLN400-HPP	Huntington's Padding Package CLN400-20, EA
------------	--