## Joerns Easytrack FS<sup>™</sup>

## Unique, Portable Ceiling Lift Solution for Residential Care

The Easytrack FS<sup>™</sup>\* is a lightweight, free-standing portable overhead track system, based on many of the same principles as the current Easytrack range. Unlike the existing Easytrack system, the Easytrack FS does not rely on support from the ceiling thereby removing any potential structural issues that may be encountered.

LIGHTWEIGHT COMPONENTS The components ensures ease of storage and transportation.

SELF LOCATING POSTS The posts are self-locating and require no additional tools or fixings.

TELESCOPIC ADJUSTMENT The Easytrack FS incorporates telescopic post and track adjustment for ease of use.

SPRING LOADED POSTS To aid set-up, the posts are spring loaded so that the height can be adjusted effortlessly.

AUTO LOCKING MECHANISM The posts automatically lock into place once adjustment has been made.

CARRY BAGS This optional accessory aid both transportation and storage.



9720010

NL	

## SAFE WORKING LOAD: 440 LB / 200 KG

SPECIFICATION		
Floor to Top of Rail	72" - 94.7"	1830-2406 mm
Floor to Bottom of Rail	66.1" - 88.8"	1680-2256 mm
Exterior Width	87" - 127.2"	2210-3230 mm
Interior Width	42.3"	1075 mm
Footplate Width	5.9"	150 mm
Rail Weight	30.9 lb	14 kg
Post Weight (each)	20.9 lb	9.5 kg
Footplate Weight	12.1 lb	5.5 kg

## **EASY TO ASSEMBLE**

The EasyTrack FS is designed to be assembled by a single person in a matter of minutes.

Specifications may vary due to manufacturing tolerances. \*Easytrack FS is a registered trademark of Arjo

