

## Features

- Scooter can be easily **disassembled into 5 pieces** for transport or storage.
- 9 mile maximum range with **standard 12AH batteries**. Range can be increased to 15 miles with **optional 21AH batteries**.
- Every unit comes standard with three sets of **replaceable color panels** allowing the user to quickly and easily switch the scooter from Metallic Red to Midnight Blue to Silver.
- Black solid, non-marking, **flat-free tires** for worry-free travel.
- Wrap-around **delta handles** make driving easy, even for those with limited dexterity.
- **Adjustable angle tiller** allows the client to position the controller into the most ideal driving position.
- **Backlit battery gauge** makes it even easier to read.
- Seat features 180 degree swivel and fold-down backrest.
- Stylish **two-tone seat upholstery**.
- Removable **flip-up armrests** have adjustable height, width and angle for a broad range of adjustments. These can easily be changed and adjusted to fit the needs of the client.
- **Armrest reflectors** make the scooter more visible when traveling in low light conditions.
- **Anti-tip wheels** for greater safety.
- Seat back folds completely flat, which comes in handy when transporting or storing the scooter.
- Standard **plastic front basket** includes removable liner bag.



# Spitfire™ EX 1420

Travel Mobility Scooter

**drive** ActiveCare™ Mobility

# Spitfire™ EX 1420 Travel Mobility Scooter

## Specifications

**HCPCS** K0800

**Weight Capacity** 300 lbs.

**Maximum Speed** 5 mph

**Maximum Range** 9 miles (12AH)  
15 miles (21AH)

**Turning Radius** 50"

**Ground Clearance** 3"

**Overall Length** 42"

**Overall Width** 20"

### Seat Dimensions

16"(w) x 15"(d) x 15"(h) Folding Seat

**Floor to Seat Height** 20"-23"

**Motor** 24V x 270W x 4700rpm

**Batteries** Two 12V x 12AH/21AH

**Battery Charger** 2 Amp, Offboard

**Brakes** Electromagnetic

**Front Wheels** 7.5" x 2" Flat-Free

**Rear Wheels** 8.5" x 3" Flat-Free

**Base Weight** 66 lbs.

**Battery Weight (pair)** 20 lbs. (12AH)  
28 lbs. (21AH)

**Seat Weight** 19 lbs.

**Heaviest Piece Wt.** 40 lbs. (Front Section)

*Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.*

## Quick & Easy Disassembly



Take your Spitfire EX Travel Scooter anywhere you want to go with its simple **disassembly feature**. Just remove the seat, front basket and battery pack, then pull up on the quick release handle located behind the seat post, and the scooter's body separates into two sections.-- perfect for loading in your car, truck or van. The sections even have built-in carrying handles!

## Inventory Numbers

SPITFIRE142016FS12.... 16" Folding Seat, 12AH Batteries

SPITFIRE142016FS21.... 16" Folding Seat, 21AH Batteries

## Accessories (sold separately)

### Carriers & Holders

- Crutch/Cane Holder
- Cylinder Tank Carrier
- Rear Basket

### Carry Bags

- Armrest Bag
- Oval Bag
- Carry All Bag

### Other Accessories

- Scooter Trailer
- Scooter Cover

## Warranty

Lifetime on frame, 14 months on electronic controller and drive train components, six months on batteries.

## Color Options

Every unit comes standard with three sets of replaceable color panels allowing the user to quickly and easily switch the scooter from Metallic Red to Midnight Blue to Silver.



**Metallic Red**



**Midnight Blue**



**Silver**

**drive** ActiveCare™ Mobility