



Deluxe Hip-High Chair

17100

Maple Wood Frame Construction, 1/cs

FEATURES:

- Designed for post-hip surgery residents.
- Provides sitting without bending at the hip.
- Armrests make sitting down or pushing up to get out from the chair safe and easy.
- Attractive vinyl seat and backrest is soft, comfortable and water resistant.
- 17100 is constructed of solid maple hardwood—not a plywood veneer.
- The fixed footrest on item # 17100 is adjustable to provide maximum comfort and an ideal ergonomic fit.
- 2 Year Warranty: Frame.
- 1 Year Limited Warranty: Upholstery.



**3 FOOTREST HEIGHT
ADJUSTMENTS!**



**BROWN VEIN
FINISH ON
17100-BV**
EASY TO
MAINTAIN & LESS
INSTITUTIONAL
LOOKING!



Hip-High Chair (Steel Construction)

17100-BV

Steel Brown Vein Frame Construction, 1/cs

FEATURES:

- Designed for post-hip surgery residents.
- Attractive, durable brown vein finish.
- Contoured armrests are comfortable and durable.
- Attractive vinyl seat and backrest is soft, comfortable and water resistant.
- 1 Year warranty on upholstery.
- Limited lifetime warranty on frame.

SPECIFICATIONS

Item #	Floor to Seat Height	Width Between Arms	Floor to Arm Height	Floor to Footrest Height (adjustable)	Overall Width	Overall Height	Overall Depth	Product Weight	Carton Shipping Weight	Weight Limit
17100	26"	19"	31"	5.5", 7.5", 9.5"	24"	41"	23.5"	33 lbs.	41 lbs.	300 lbs.
17100-BV	26"	17.5"	32"	6"	22"	40.75"	29"	33.6 lbs.	41.6 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

patient room

drive
WWW.DRIVEMEDICAL.COM

To order call toll free: 877.224.0946