

Premier LifeCare Recliner 5400



ISO 9001-2008 Certified

The Winco 5400 is an upscale long-term recliner with fully upholstered sides for contemporary modern style, and practical function. Durable and comfortably contoured urethane armrests protect the upholstered finish and are easily cleaned. The 5400 is a perfect blend of value, function and style.



KEY FEATURES

- Urethane armrests
- Fully upholstered
- Multiple side-table Options
- 3-position lockbar
- "No pinch" safety zone

Proudly Made in the USA




1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

PREMIER LIFE CARE RECLINER 5400



- 1 Push-handle for easy maneuverability.
- 2 Comfortable, contoured urethane armrests to keep upholstery free from damage.
- 3 Low-profile 3" locking casters provide smooth & easy maneuvering.
- 4 3 position lock-bar feature allows for safe, secure, and comfortable positioning.



SPECIFICATIONS

5400


Seat Height	20.5" (52cm)
Seat Width	19.5" (49.5cm)
Seat Depth	20" (50.8cm)
Overall Height	46" (116.8cm)
Overall Width	27.75" (70.5cm)
Back Height Above Seat	30.5" (77.5cm)
Overall Length (Upright)	39.5" (100.3cm)
Overall Length (Reclined)	62" (157.5cm)
Weight	90 lbs. (40.8kg)
Weight Capacity	350 lbs. (158.7kg)

Add 1" to chair width for each side-table in down position
 *Add 10.5" to chair width for each raised side table

OPTIONS

IB-133 TB	IV Pole IV	Oxygen Tank Holder O2	Pivot-Table PN	Designer Black Wheels (Note: Lowers Arm, Seat and Overall Heights by 0.75") Model 5405
Side-Table LT, RT, DT	Loriso Support T500	Contoured Headrest CH00	Lateral Support LS00	Side Cushion SC00


STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

OPTIONAL COLOR UPGRADES

(ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL

(ADDITIONAL CHARGES APPLY)



COLORS

10-03-11

REV.D

009100