

LifeCare Recliner 5851/585S/5861/586S



ISO 9001-2008 Certified

Value, durability, and comfort in an exclusive recliner design

Unsurpassed quality, durability, and design are standard for Winco's basic recliner. Blow-molded side panels, built-in hand-grips, total-lock steel casters, and round, powder-coated steel frame make the LifeCare Recliner a perfect choice for value.

The 585S/586S includes a highly specialized "saddle" seat, designed to prevent skin breakdown for patients. A polyurethane seat cover provides better airflow while an internal gel-pack molds to the body for extra protection.

KEY FEATURES

- Non-porous tray
- Plastic molded side panels
- Padded armrest with arm guard
- "No pinch" safety zone



5851 - w/ tray

5861 - w/out tray

US Patent #7070243



Available with "saddle" seat

585S - w/ tray

586S - w/out tray

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474


www.wincomfg.com



PRODUCT SPECIFICATIONS

LIFECARE RECLINER 5851/585S/5861/586S

- 1 Padded push-handle for easy maneuverability.
- 2 Comfortable, padded armrests with guards to keep upholstery free from damage.
- 3 Non-porous tray, stores on side of chair.
- 4 5" Total lock casters (in rear) allow for better handling and control.
- 5 3 position lock-bar allows for safe, comfortable positioning.
- 6 Plastic side panels help to create a "safe zone" for patient care without pinch points.
- 7 585S/586S "Saddle" seat model provides extra care and helps to prevent skin breakdown.



SPECIFICATIONS

| | 5851/5861 | 585S/586S |
|---------------------------|-------------------------------|-------------------------------|
| Seat Height | 20" (50.8cm) | 19" (48.3cm) |
| Seat Width | 19" (48.3cm) | 19.5" (49.5cm) |
| Seat Depth | 20" (50.8cm) | 21" (53.3cm) |
| Overall Height | 46.5" (118.1cm) | 46.5" (118.1cm) |
| Overall Width | 26.5" (67.3cm) | 26.5" (67.3cm) |
| Back Height Above Seat | 32" (81.3cm) | 31.5" (80cm) |
| Overall Length (Upright) | 41" (104.1cm) | 41" (104.1cm) |
| Overall Length (Reclined) | 63" (160cm) | 63" (160cm) |
| Weight | 5851: 85 lbs. (38.6kg) | 585S: 85 lbs. (38.6kg) |
| | 5861: 80 lbs. (36.3kg) | 586S: 80 lbs. (36.3kg) |
| Weight Capacity | 275 lbs. (124.7kg) | 275 lbs. (124.7kg) |

OPTIONS



TB-133 TB



LiquiCell LC
(5851/5861 ONLY)



IV Pole IV



Flip-up Footrests
F500



Torso Support
T500



Contoured Headrest
CH00



Lateral Support
L500



Side Cushion
SC00

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



COLORS