

Med - Aire Plus Bariatric Alternating Pressure Mattress Replacement System

14048	48" (W) x 80" (L) x 10" (H), Pump and Mattress System, 1/cs
14048M	48" (W) x 80" (L) x 10" (H), Mattress Only for 14048
14030XP	Pump Only for 14048
14048-CELL	Replacement Bladder for 14048, 1/ea

HCPCS code: E0277

FEATURES:

- The product is designed to be used in the prevention, treatment and management of Pressure Ulcers.
- Mattress operates in an alternating pressure mode.
- Ideal for use in nursing homes, hospitals and home use.
- LED and Audible low pressure alarm activate when pressure falls below 10-mmHg.
- Cell-on-Cell design prevents the patient from "Bottoming Out" in the event of power failure.
- "Whisper quiet" pump produces 12-LPM (Litres per Minute) in air flow.
- Mattress has 20 air bladders.
- **10" cells are constructed with 60% Nylon/40% PU that is soft, lightweight and will not stretch.**
- Quick fill adapter standard.
- CPR valves provide fast (20 second) deflation.
- The zippered Nylon/PU cover material is quilted, water resistant, vapor permeable and machine washable.

SPECIFICATIONS:

Pump:

- Power: 120 VAC
- Input Frequency: 60 Hz
- Pressure Range: 20 mmHg to 60 mmHg (inner mattress pressure)
- Auto firm (60 mmHg, Auto return in 30 min.)
- Seat function (Pressure preset + 10 mmHg)
- Static Selectable
- Static auto return (in 2 hours)
- Comfort therapy
- Function lock (Auto lock in 30 min. after on)
- Alarm mute
- Low pressure alarm
- System (alternating) passive alarm
- Power failure alarm
- Low pressure alarm delay
- Memory function
- Air Flow: 12-LPM (Litres Per Minute)
- Functions: Static and dynamic features
- Power Cord Length: 15-Ft.
- Pump Dimensions: 14.2" (W) x 9.8" (H) x 4.5" (L)
- UL: 60601-1, 1st Edition-2003
- Maximum Current: 0.5 Amps
- Warranty: 18-Months Limited

Mattress Dimension:

- Topper Material: Zippered and Quilted Urethane Coated Nylon Fabric (2-Way Stretch) is; water resistant, vapor permeable, machine washable and low shear
- Base Material: Zippered Vinyl
- 4-Air cells at head section remain static to enhance pillow function
- Dimensions: 48" (W) x 80" (L) x 10" (H)
- Maximum Support Capacity: 750 lbs.
- Warranty: 18-Months Limited



BARIATRIC
750 LB
WEIGHT CAPACITY

PART # 14427A NOW COMES WITH "QUICK FILL VALVE" FOR USE WITH PART # 14427, QUICK INFLATE ELECTRIC PUMP (SOLD SEPARATELY)



14427 (Electric Pump)

**NEW
DIGITAL
PUMP
DESIGN!**



14030XP
**PUMP ONLY
NEW & IMPROVED!**



14007

EZ Lock/Quick Release Tubing Connector

- Standard with 14025, 14027, 14028N, 14029, 14029DP, 14030, 14048, 14054 and 14060.
- Locks pump and mattress together.
- Pressure tabs allow easy and quick disconnect.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.