



AEROSOL DELIVERY SYSTEM



VIOS® AEROSOL DELIVERY SYSTEM

- ✓ Highly efficient & powerful compressor system
- ✓ FAST treatments – 5-6 minutes¹
- ✓ Saves money by wasting less medication
- ✓ Delivers consistent particle size with targeted delivery to the lungs
- ✓ Includes 2 PARI Reusable Nebulizers
 - Easy to clean – hand wash, boilable and dishwasher safe
 - Reusable nebulizer is designed to last 6 months
 - Covered by insurance every 6 months²

PARI. The Clinical Choice.

PARI VIOS®

INCLUDING PARI TIMESTRIP®



Activated


6 Months Later

- ✓ Easy to use visual indicator
- ✓ Reminds patients to replace their PARI nebulizer every 6 months
- ✓ Helps maintain FAST treatments & reduce bacterial contamination³

THE DIFFERENCE IS IN THE DELIVERY

**PARI Reusable Nebulizers
vs.
Disposable Nebulizers**

SLOW ← ————— → **FAST**



Faster treatments with PARI⁵

CONFIGURATIONS

	Part No.	NDC No.	HCPCS No.
Standard, LC® Sprint	310F35-LCS	83490-0310-58	E0570, A7005
Standard, LC® Plus	310F83-LC+	83490-0310-67	E0570, A7005
Deluxe, LC® Sprint	310F35-LSD	83490-0310-59	E0570, A7005
Deluxe, LC® Plus	310F81-LCD	83490-0310-64	E0570, A7005

SPECIFICATIONS

Technical Specifications	
Maximum Pressure	34 PSI (2.3 bar)
Maximum Liter Flow	8.5 LPM
Nominal Operating Pressure ⁴	18 PSI (1.2 bar)
Nominal Air Flow ⁴	4.5 LPM
Weight	3 lbs (1300 gms)
Dimensions (Length, Width, Height)	6.5 x 6.5 x 3.75 in
Compressor Warranty	5 years (6 years user activated)

Accessories	
PARI LC® Sprint	023F35
PARI LC® Plus	022F81
PARI Adult Aerosol Mask	044F7251
Replacement Air Filters	041E4851P12



All configurations include a PARI Reusable Nebulizer and backup reusable nebulizer, Vios® Compressor, Timestrip®, and Wing-Tip™ tubing. Deluxe includes Instructional DVD and Dual Compartment Carrying Case.

¹Measured with Malvern MasterSizer X at 50% relative humidity, 0.9% NaCl solution, inspiratory flow 20 liters/minute, continuous nebulization, 23° C, fill volume 2.5 ml.

²Insurance coverage varies by provider

³Lester, et al. Nebulizer Use and Maintenance by Cystic Fibrosis Patients: A Survey Study. *Respiratory Care* 2004; 49(12): 1504-1508.

⁴When used with a PARI reusable nebulizer.

⁵Tiffin, et al. Comparison of Nebulized Drug Delivery Times of Four Nebulizers. *CSRT* 2007.