



AEROSOL DELIVERY SYSTEM



VIOS[®] AEROSOL DELIVERY SYSTEM

- ✓ Highly efficient & powerful compressor system
- ✓ FAST treatments – 5-6 minutes¹
- ✓ Saves money by wasting less medication
- ✓ Delivers consistent particle size with targeted delivery to the lungs
- ✓ PARI Pals™ pediatric extras makes treatments fun for kids!
- ✓ Includes 2 PARI Reusable Nebulizers
 - Easy to clean – hand wash, boilable and dishwasher safe
 - Reusable nebulizer is designed to last 6 months
 - Covered by insurance every 6 months²

PARI. The Clinical Choice.

PARI VIOS®

INCLUDING PARI TIMESTRIP®



Activated

6 Months Later

- ✓ Easy to use visual indicator
- ✓ Reminds patients to replace their PARI nebulizer every 6 months
- ✓ Helps maintain FAST treatments & reduce bacterial contamination³

PARI PALS™ MAKE TREATMENTS FUN FOR KIDS!

PARI PALS
PLAYGROUND



Visit PARIPals.com

CONFIGURATIONS

	Part No.	NDC No.	HCPCS No.
Pediatric, LC® Sprint	310F35-P	83490-0310-62	E0570, A7005, A7015
Pediatric, LC® Plus	310F83-P	83490-0310-71	E0570, A7005, A7015

SPECIFICATIONS

Technical Specifications

Maximum Pressure	34 PSI (2.3 bar)
Maximum Liter Flow	8.5 LPM
Nominal Operating Pressure ⁴	18 PSI (1.2 bar)
Nominal Air Flow ⁴	4.5 LPM
Weight	3 lbs (1300 gms)
Dimensions (Length, Width, Height)	6.5 x 6.5 x 3.75 in
Compressor Warranty	5 years (6 years user activated)

Accessories

PARI LC® Sprint	023F35
PARI LC® Plus	022F81
Bubbles the Fish™ II pediatric aerosol mask	044F7248
Replacement Air Filters	041E4851P12



All configurations include a PARI Reusable Nebulizer and backup reusable nebulizer, Vios® Pediatric Compressor, Timestrip®, Bubbles the Fish™ II pediatric aerosol mask, Wing-Tip™ tubing, PARI Pals™ Activity Book & Sticker Set, PARI Pals™ Decoder Card, Instructional DVD and Dual Compartment Carrying Case.

¹Measured with Malvern MasterSizer X at 50% relative humidity, 0.9% NaCl solution, inspiratory flow 20 liters/minute, continuous nebulization, 23° C, fill volume 2.5 ml.

²Insurance coverage varies by provider

³Lester, Mary K. Nebulizer Use and Maintenance by Cystic Fibrosis Patients: A Survey Study. *Respiratory Care* 2004; 49(12): 1504-1508.

⁴When used with a PARI reusable nebulizer.