

## SPECIFICATIONS for Drive Nimbo Rehab Lightweight Posterior Posture Walker

\*FOR PROPER SIZING: Measure wrist to floor height

Item #	Handlebar Height Range	Overall Width	Width Inside Handle Grip	Depth Opened at Base	Product Weight	Weight Capacity
KA 1200N (Tyke)	14.5" - 18"	24"	13.5"	25"	9.5 lbs.	100 lbs.
KA 2200N (Youth)	19" - 25"	24"	13.5"	31"	10.5 lbs.	100 lbs.
KA 3200N (Junior)	23" - 30.5"	26"	15"	34.5"	12.5 lbs.	200 lbs.
KA 4200N (Yg. Adult)	28"-36"	26.5"	16.5"	37"	14 lbs.	200 lbs.