



**TRANSPORT CHAIR &
ROLLATOR IN 1 PRODUCT!**



**COMES WITH
PADDED SEAT
& POUCH!**

**FOLD DOWN FOOT REST FOR
USE AS TRANSPORT CHAIR!**



Duet Transport Chair/Rollator

795B	Blue, 1/cs
795BK	Black, 1/cs
795G	Gold, 1/cs
795BU	Burgundy, 1/cs

FEATURES:

- Combines the features of a transport chair and rollator so an individual can ambulate by themselves or be pushed by a caregiver.
- The two position padded contoured back rest can be attached to the front or the back provides maximum comfort and convenience.
- Weighs only 19 lbs.
- Large 8" casters are ideal for indoor and outdoor use.
- Easy release loop locks.
- Comes with carry pouch.
- Seat lock prevents the Duet from accidentally folding
- The Duet comes with two sets of spacers to raise or lower the seat height.
- The foot rest folds down when used as a transport chair and folds up and out of the way when used as a rollator.
- Comfortable padded seat.
- Limited Lifetime Warranty.
- **Lifetime Warranty on brake cable.**

**NEW LOCK LOCATED ON
BOTTOM OF SEAT!**



**NEW SEAT LOCK
PREVENTS ROLLATOR
FROM ACCIDENTALLY
FOLDING.**



**EASY HEIGHT
ADJUSTMENT
WITH SELF
THREADING
KNOB.**

SPECIFICATIONS

Item #	Overall Width (open)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length w/Riggings	Product Weight	Carton Shipping Weight	Weight Limit
DUET	23"	36"	24" - 26"	36"	12"	12"	35"	19 lbs.	25 lbs.	300 lbs.

drive MEDICAL DESIGN & MANUFACTURING

rollators

To order call toll free: 877.224.0946