



# Aqua Creek Products

The Total Solution for Pool and Spa Access

## The Portable Pro Pool™ Lift

U.S. Patent No.  
5790995

Item# F-004PPPB

The Portable Pro Pool™ Lift is a completely portable aquatic lift that can accommodate loads up to 400lb (181kg). It requires no anchors and can easily be removed from poolside when not needed.

The lift features a compact, low profile frame that does not intrude into the pool when not being used. It's aesthetically pleasing and compact design make it a popular lift.

### Features:

- 400 Pound (181 kg) weight capacity
- Rechargeable battery with wall mount charger
- Adjustable footrest & lap belt
- Dual flip-up arm rests
- Simple push button operation
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty

### Options:

- Upgrade Pack- Item #F-004PPPBAP  
Includes: battery, headrest, chest strap and cover
- Headrest - Item #F-422PLH
- Chest strap - Item #F-423CS
- Protective cover - Item #F-444PPPC
- Cycle Attachment - Item #F-019CA



**PROUDLY SOLD BY:**