



FEATURES

- Durable, low-maintenance, triple chrome-plated, carbon steel frame is long-lasting
- Dynamic recline range from 90° to 168°
- Rear axle base moves, providing stable center of gravity in recline position
- Multiple frame color variations for customizing
- Multiple wheel offerings allow variation of seat-to-floor heights
- Urethane rear tires mounted on “no flex” wheels offer superior performance
- Urethane casters for better ride and lighter weight
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- Additional “Lift Packages” assist attendants in positioning patients
- Stocked with attendant (standard) or patient-operated wheel locks
- Meets RESNA standards for manual wheelchairs in product performance

Invacare® IVC™ 9000 XT Recliner

The Invacare IVC 9000 XT Recliner wheelchair is the number one choice for individuals who demand a truly versatile, yet moderately priced lightweight wheelchair. The 9000 XT Recliner can be configured into an endless number of seat widths, depths and heights, and has multiple frame and upholstery options. The recliner has been standardized to interchange components and accessories with the other 9000 series chairs. Whether you need the lighter recliner package or the assist recliner package, the 9000 XT Recliner can provide it.

Specifications

Dimensions

Seat sizes

Widths	14", 16", 18", 19", 20", 22"
Depths	16", 17", 18"
Seat-to-floor heights	19.5" 17.5", 18.5" with rear wheel package
Back height	32" recliner (22" high back, 10" removable headrest)

Front riggings

Hemi swingaway footrest	17"-20"
Hemi swingaway padded elevating legrest	17.5"-21"
Hemi articulating elevating legrest	17"-25"

(Footplates available in aluminum or composite)

Arm styles

Full-length fixed height, padded, removable	
Desk-length fixed height, padded, removable	
Full-length adjustable height, padded, removable	
Desk-length adjustable height, padded, removable	

(All arms are conventional)

Weight capacity	250 lb.
Product weight (9RC/WD86/R/28/2R)	51 lb. without lift package 62 lb. with lift package
Warranty	lifetime on frame and crossbraces