



Invacare® IVC™ Tracer® SX5 Recliner

The Invacare IVC Tracer SX5 Recliner wheelchair offers the durability of a lightweight manual wheelchair frame with the added versatility and comfort of a full reclining back. The Tracer SX5 Recliner offers dynamic recline range from 100° to 168°. The Tracer SX5 Recliner offers a 300 lb. weight capacity on 20" and 22" widths. This chair also is available with multiple seat widths, depths, armstyles, front riggings, a full range of accessories and a lifetime warranty.

FEATURES

- The Tracer SX5 Full Recliner offers a dynamic range of recline from 100° to 168°
- Durable, low-maintenance, triple chrome-plated, carbon steel frame
- Urethane rear tires, mounted on "no flex" wheels, offer superior performance
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- Upholstery is dual-embossed
- Meets RESNA standards for manual wheelchairs in product performance

Specifications

Dimensions	
Seat sizes	
Widths	14", 16", 18", 20", 22"
Depths	16" standard, 18" by seat extension (optional)
Seat-to-floor height	19.5"
	17.5", 18.5" with different wheel and caster packages
Back height	33" recliner (23" high back, 10" removable headrest)
Front riggings	
Hemi swingaway footrest	17"-20"
Hemi swingaway padded elevating legrest	17.5"-21"
Hemi articulating elevating legrest	17"-25"
(Footplates available in aluminum or composite)	
Arm styles	
Full-length fixed height, padded, removable	
Desk-length fixed height, padded, removable	
Full-length adjustable height, padded, removable	
Desk-length adjustable height, padded, removable	
(All arms are conventional)	
Weight capacity	250 lb. (300 lb. 20" and 22" width)
Product weight (TRSX5/WD86/RC/28)	53 lb.
Warranty	lifetime on frame and crossbraces