



Jazzy[®]

ELITE 14

The Jazzy[®] Elite 14 offers a 300 lbs. weight capacity and excellent performance in a stylish, highly maneuverable package.



Pride[®]
Mobility Products Corp.
LIVE YOUR BEST[®]

SPECIFICATIONS

Weight capacity	300 lbs.
Maximum speed	Up to 4 mph
Ground clearance	1.5" at anti-tips
Turning radius	24" without front riggings
Overall length	35.5" without front riggings
Overall width	24.25"
Drive wheels	14"
Front wheels	3" solid anti-tips
Rear wheels	6" solid casters
Drivetrain	Two-motor, front-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	50 amp controller
Specialty controls	N/A
Battery charger	Off-board, 5 amp (standard)
Per-charge range	9.8 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	NF-22
Battery weight	37 lbs. each
Base weight	98.5 lbs.
Standard seat weight	51.5 lbs.
Max. seat dimensions	22" x 22"
Warranty	5-year limited warranty on frame; 13-month limited warranty on electronics; 13-month limited warranty on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

MODELS

Jazzy Elite 14 2S-SS	(Code: K0822)
Jazzy Elite 14 2S-C	(Code: K0823)

FEATURES

- Two-motor, front-wheel drive
- 14" knobby tires and optimum weight distribution provide superior traction
- NF-22 batteries for maximum range
- 5 amp, off-board charger

OPTIONS

- High-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

RECOMMENDED LIFTS

Pride® Lifts & Ramps

- Backpacker® Plus
- Outlander
- Backpacker® MV



COLOR CHOICES



Red



Blue

PRIDE[®]
PERFORMANCE MOBILITY
— LIVE YOUR BEST[®] —

(US) 800-800-8586 ▪ (Canada) 888-570-1113 ▪ www.pridemobility.com