




Kanga TS Pediatric Tilt-In-Space Wheelchair

Available
1 Pack

KG 1000 Kanga TS 10", Yellow, 1/bx

KG 1200 Kanga TS 12", Red, 1/bx

KG 1400 Kanga TS 14", Blue, 1/bx



NEW!

AVAILABLE
IN 10", 12"
AND 14"
SEAT
WIDTHS.



FEATURES:

- Aluminum frame.
- Features an economical seating system that includes depth adjustable seat, five-point harness, depth adjustable abductor, hip guides, lateral supports and curved headrest.
- Height and angle adjustable flip-back armrests.
- Height adjustable push handles.
- Tilt-in-Space up to 35 degrees.
- Seat to back angle adjustability.
- Angled seat to prevent forward sliding.
- Swing away, height adjustable plastic footrests.
- Angle adjustable plastic footplates.
- 20" rear mag wheels with semi-pneumatic quick release wheels.
- Frame adapts to most custom rehab seating and positioning systems.
- Crash tested for bus transit.
- Folds for easy transport.
- Weight Capacity 200 lbs.
- Limited Lifetime Warranty.



Red



Blue



Yellow

SPECIFICATIONS

Seat Width	Overall Width (open)	Overall Width (closed)	Front Seat To-Floor Height	Back Seat To-Floor Height	Seat Depth	Overall Height w/ Headrest	Overall Depth w/ Riggings	Product Weight	Carton Shipping Weight	Weight Limit
10"	18"	14"	21"	20"	13"-15"	38.9"-41"	37"	36.5 lbs.	48.5 lbs.	200 lbs.
12"	20"	14"	21"	20"	13"-15"	38.9"-41"	37"	36.5 lbs.	48.5 lbs.	200 lbs.
14"	22"	14"	21"	20"	13"-15"	38.9"-41"	37"	36.5 lbs.	48.5 lbs.	200 lbs.

mobility
wheelchairs/kanga ts

WenzeliteRe/hab

To order call toll free: 877.224.0946

drive