

Luggie Scooter Models

FEATURES	LUGGIE (STANDARD)	LUGGIE (DELUXE)	LUGGIE (ELITE)
Seat Size – width (inches)	Standard (14.5")	Larger (16.5")	Larger (16.5")
Set of Armrest	Optional	Optional	Included
Adjustable Steering Heights	Yes	Yes	Yes
Speed Control Dial	Yes	Yes	Yes
Battery Indicator	Meter	Meter	LED
Retractable Handlebar	No	No	Yes
Auto Retractable Bumper	No	No	Yes
2 Seat Height Adjustments (3 locking positions)	Yes	Yes	Yes
Overall Length (unfolded)	39"	39"	39"
Overall Length (folded)	26"	26"	26"
Overall Width	18"	18"	18"
Ground Clearance	2.5"	2.5"	2.5"
Weight Capacity	250 lbs	250 lbs	320 lbs
Range*	11 mi	11 mi	6.87 mi**
Total Weight without Battery	51.7 lbs	51.7 lbs	51.7 lbs
Turning Radius	35.5"	35.5"	35.5"
Maximum Speed	4 MPH	4 MPH	4 MPH
Lithium Battery Strength	8.5 AH	10.5 AH	10.5 AH
Alloy aluminum frame	Yes	Yes	Yes
Rear wheel drive	Yes	Yes	Yes
24 Volt DC magnetic motor	Yes	Yes	Yes
Intelligent braking system	Yes	Yes	Yes

*Range may vary with weight, terrain, temperature, battery condition, etc. and is subject to manufacturing tolerances. Charge your battery when the indicator is in the red zone on your battery indicator instead of letting the power run out, for longer life of your battery.

**Range obtained from testing on (146 kg) 321 pound capacity. Range may vary with weight, terrain, temperature, battery condition, etc. and is subject to manufacturing tolerances. Charge your battery when the indicator is in the red zone on your battery indicator instead of letting the power run out, for longer life of your battery.

