

zip'r Mantis



**Seat
Elevator
Option
Available**

The Zip'r Mantis power wheelchair features a true mid-wheel design that maximizes maneuverability, while caster wheels in both the front and rear offer premium stability. Electromagnetic, regenerative brakes, the highly programmable PG VR-2 controller, and flat-free tires are some of the many standard features included with the Zip'r Mantis. Upgrade to a Zip'r Mantis with a seat elevator for five inches of seat-height control with the touch of a joystick.

- 6" front and rear casters offer premium stability
- Mid-wheel design maximizes maneuverability
- Flat-free tires
- Up to 15 mile range per battery charge
- Highly programmable, 50 amp, PG VR-2 controller
- 300 lb. weight capacity

Specifications

Weight capacity	300 lbs.
Maximum speed*	4.9 mph
Per-charge range*	Up to 15 miles
Turning radius	22.4"
Ground clearance	2"
Overall length	42.3"
Overall width	26"
Overall height	49" - 51" 49" - 54" (SE)
Weight w/ batteries	189 lbs.
Weight: Base	101 lbs.
Weight: Seat	43 lbs.
Weight: Batteries	45 lbs.
Heaviest piece	101 lbs. (the base)
Battery	2 U1 (12 V x 35 Ah)
Charger	Off-board, 5A
Drive	Mid-wheel
Motor specifications	200 W (.27 hp), 24 V
Controller	PGDT VR-2, 50A
Brakes	Electromagnetic, regenerative
Seat type	Captain seat
Seat width	18"
Seat depth	18"
Seat back height	30"
Seat-to-floor height	20" - 22" 21.7" - 26.8" (SE)
Drive wheel	10" x 3.5" flat-free
Front caster wheel	6" x 2" flat-free
Rear caster wheel	6" x 2" flat-free
Anti-tip wheel	N/A
Body colors	Red/Blue
HCPCS code	K0823

Limited Warranty**

Frame	3 years
Drivetrain	1 year
Electronic Components	1 year
Batteries	6 months

All specifications are subject to change without notification.

*Specification may vary due to user weight, speed, battery condition and terrain.

**Exclusions apply. See owner's manual for more details.