## **Metro Features**

- 30° Fixed Tilt Increases upper body positioning for someone with limited upper body strength
- Lightweight & Compact Folding Frame Offers unbeatable portability – easy to stow and transport
- New Upholstery
   Has added contouring to provide greater
   support and comfort
- Mesh Upholstery Center Increases ventilation to keep the occupant cooler
- Under Seat Storage Basket Is the perfect option for storing essential items

- Adjustable Trunk Support Promotes superior lumbar positioning
- New Wider Wheels Provide added stability and a comfortable ride on uneven terrain
- "Quick Release" Rear Wheels Easily remove for storage or transport
- Rigorous Product Testing Ensures maximum safety, strength, stability and endurance
- Crash Tested Transit Models Adhere to national safety standards, and provide safe transport in bus or van



Positioning options and accessories sold separately.

## **Metro Specifications**

Measurement	Chair Model	Seat Width*	Seat Depth*	Seat Back Height	Seat to Footplate*	Seat to Floor*	Seat Angle	teadrest Extension	Weight Capacity Regular/Transit	Weight of Chair <sup>A</sup>
in.	ME 12	12	11	21	8 - 15	24	30°	6	75 lbs / 66 lbs	23 lbs
cm	ME 12	30.5	28	53.3	20.3 - 38	61		15	34 kg / 30 kg	10.4 kg
in.	<b>ME</b> 14	14	13	23	8 - 14.5	24	30°	8	100 lbs / 100 lbs	24 lbs
cm	<b>ME</b> 14	35.5	33	58.4	20.3 - 36.8	61		20	45.5 kg / 45.5 kg	11kg
in.	ME 16	16	14	25	13 - 18	23	30°	8	170 lbs / 170 lbs	27 lbs
cm	ME 16	40.5	35.6	63.5	33 - 45.7	58.4		20	77 kg / 77 kg	12.3 kg
in.	ME 18	18	16	25	3 - 24	23.5	30°	8	250 lbs / 170 lbs	29 lbs
cm	ME 18	45.5	40.6	63.5	7.6 - 61	59.7		20	114 kg / 77 kg	13.2 kg

in. - Inches, cm - Centimeters \* Plus or minus .5 inch (1.27 cm) due to fabric. Δ - Weight of chair is without front rigging.

Measurement	Chair Model	Overall Width	Overall Height	Overall Length	Folded Length	Folded	Folded Height
in.	ME 12	23.5	36.75	36.5	48	15	14
cm	ME 12	60	93.3	92.7	122	38.1	35.6
in.	<b>ME</b> 14	25.5	37.5	34.5	48	15.5	14
cm	<b>ME</b> 14	64.8	95.3	87.6	122	39.4	35.6
in.	ME 16	27.5	38	35	49	17	14.5
cm	ME 16	70	96.5	89	124.5	43.2	37
in.	ME 18	28.5	38	42	51	14.75	13.75
cm	ME 18	72.4	96.5	106.7	129.5	37.5	35

Use only Convaid options and parts on Convaid products. Convaid parts are not interchangeable with other manufacturer's products. Replace any worn parts immediately.

Note: Convaid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.


## **Limited Lifetime Warranty**

Frame and X-Braces......Lifetime Fabric and webbing......One year Other components.....One year

This warranty does not cover normal wear and tear, or damage caused by accident or misuse.



## **Upholstery Colors**



CE

U.S. Patent 3,995,882 This product may be covered by the following U.S. patent - 6,099,075 Additional patents pending.

Rev 6-29-2010

visit us at www.convaid.com