

H E A V Y D U T Y P O V



Seat Color
Charcoal Vinyl

Body Colors
Orange
Red

Ride

Golden Patriot™ 4 wheel

If you're looking for the ultimate outdoor driving experience, look no further than the Golden Patriot™. The Golden Patriot's rugged outdoor design features an off road package, four-wheel independent suspension, a full lighting package and a full sized captain's seat. With a 400 lb. weight capacity and 13" solid, low profile tires, the Golden Patriot can give you the ride of your dreams.



Rise. Roll. Ride.



YOUR LIFE IN MOTION

Patriot 4 wheel Specifications

Model	GR575
Medicare Code	Not Coded
Weight Capacity	400 lbs
Drive Wheels	Rear
Maximum Speed	8 mph
Operating Range ¹	13 miles
at 200 lbs	18 miles
Ground Clearance	3.5"
Turning Radius	76.5"
Type Batteries	2-75 AH
Colors	Red, Orange
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Mechanical Hand Brake	Yes
Charger	Off Board
Controller	PG S-140A
Dimensions:	
Length	57"
Width	28"
Height (ground to top of back of seat)	46.5" - 48.5"
Height (ground to top of back of headrest)	54" - 59"
Ground to Top of Seat	27" - 29"
Ground to Top of Deck	8.75"
Deck to Top of Seat	18" - 20"
Number of Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	40"
Weight of Unit: (assembled)	314 lbs
Base	160 lbs
Seat w/Arms	50 lbs
Batteries	52 lbs ea
Tires:	Solid, low profile
Front	13"
Rear	13"
Standard Seat:	Captain's Highback
Slide Seat	N/A
Back Height - no headrest	21"
Back Height with headrest	28.5" - 31.5"
Width x Depth (adjustable depth is standard, set to 18")	20" x 18"/20"
Color	Charcoal Gray
Warranty:	
Frame	5 Years
Drive Train	1 Year
Electronics ²	1 Year

¹ Battery range will vary due to rider weight, drive surface and terrain.

² Electronics warranty excludes batteries.



Full front and rear lighting package standard




Floor deck features rugged rubber mat and anti-skid material



Four-wheel independent suspension standard



Water bottle holders for cold beverage bottles only. Removable for easy cleaning.



Control panel features safety horn, flashers, lights, turn signals, speed knob and backlit battery gauge



Black, low profile, non-marking solid tires with chrome hubcaps standard



Hand brake standard