

**Neb-u-Tyke Speedster**

- Colorful, Child-Friendly Racecar design
- Interactive stickers for patient to customize their car
- Reliable Piston-driven Compressor
- Comes complete with Nebulizer Accessory Kit

<b>JB0112-067</b>	Neb-u-Tyke Speedster	1 ea
<b>JB77200</b>	Neb-u-Tyke Speedster with Garage	1 ea
<b>JB76220</b>	Garage Bag only	1 ea

**The John Bunn Neb-u-Tyke Speedster comes complete with:**

(single patient use only)

- 1 reusable nebulizer
- 1 adult mask
- 1 pediatric mask
- 1 mouthpiece
- 1 PVC tubing
- 1 "T" adapter
- 3 spare filters



**Technical Specifications**

Power .....	115VAC, 60Hz
Power Consumption .....	≤ 65 watts
Operating Temperature Range .....	41°F TO 104°F
Operating Humidity Range .....	30 to 95% RH
Storage Temperature Range .....	14°F TO 158°F
Storage Humidity Range .....	10 to 95% RH
Width .....	6.9"
Height .....	3.7"
Depth .....	10.2"
Weight .....	3.8 lbs
Max Pressure (Compressor) .....	30 psi
Flow Range (Compressor) .....	0-12 liters/min
Particle Size (Nebulizer) .....	≤ 5 microns
Capacity (Nebulizer) .....	5 ml
Method Of Nebulization .....	Piston Driven



**Graham-Field Health Products**

2935 Northeast Parkway  
Atlanta, Georgia 30360

2007© Graham-Field Health Products

For Placing Orders:

**1-800-347-5678**

**Fax: 1-800-726-0601**

**www.grahamfield.com**

GF080003RevA08