

REBEL

NEW!



CARTON SIZE MAKES THE REBEL
PERFECT FOR RETAIL!



PACKAGED IN SMALL CARTON THAT
CAN BE SHIPPED UPS!

Rebel Wheelchair

Single Axle

FEATURES:

- Quick release rear wheels allow for easy storage or transport.
- Fold-down back and folding frame allows transport in almost any trunk.
- Easily removable swing away footrests are also convenient for storage.
- Flip up desk style armrests make transfers easy and accessibility better.
- Lightweight and extremely maneuverable, easy to navigate tight areas.
- Large accommodating seating surface.
- Limited Lifetime Warranty.

QUICK
RELEASE
REAR AXLE



FOLD
DOWN
BACK



TOOL FREE
REMOVABLE
FRONT
CASTERS



1-Pack Available for All Wheelchair Models

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

Rebel – SINGLE AXLE

RTLREB18DDA-SF	18", Flip Back Desk Arm, Swing-away Footrest	1/cs
----------------	--	------

REPLACEMENT FRONT RIGGINGS

RTLREB-12	Swing-Away Footrest	1pr/cs
-----------	---------------------	--------

ARMS

RTLREB-20	Right, Flip Back Desk Arm	1/cs
RTLREB-19	Left, Flip Back Desk Arm	1/cs

WHEEL LOCKS (Push to Lock)

STDS8M241R	Right Wheel Lock, Push to Lock, Detachable Arm	1/ea
STDS8M241L	Left Wheel Lock, Push to Lock, Detachable Arm	1/ea

SEAT BELTS

STDS850 (Auto Clasp Type)	1/bx
STDS851 (Velcro® Type Closure)	1/bx

BARIATRIC SEAT BELTS

STDS855 (Auto Clasp Type)	1/bx
STDS856 (Velcro® Type Closure)	1/bx

UNIVERSAL OXYGEN "E" CYLINDER HOLDER

STDS803	
STDS804 with I.V. Pole Attachment	

WHEEL LOCK EXTENSIONS

STDS801	1/ea
---------	------

ANTI FOLD BAR

STDS806	1/cs
---------	------

WHEELCHAIR CANE/CRUTCH HOLDER

STDS1034	1/bx
----------	------

HEEL LOOPS

STDS831	1pr/bx
---------	--------

ACCESSORIES ON PAGES 57-78**SPECIFICATIONS****1-Pack Available for All Wheelchair Models**

Seat Width	Overall Width (open)	Overall Width (closed)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length (with riggings)	Product Weight	Carton Shipping Weight	Weight Limit
18"	25.5"	10.5"	35"	19.5"	17"	18"	16"	43"	37 lbs.	42 lbs.	250 lbs.

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.