

Aqua Creek Products

The Total Solution for Pool and Spa Access



The Revolution™ Lift




Gurney Option



Spineboard Option

OPTIONS

- **Upgrade Pack - Item #F-700RLAP**
Includes battery, headrest, chest strap, cover
- **Anchor, standard** (additional) -
Item #F-710RLA
- **Headrest - Item #F-722HRR**
- **Chest strap - Item #F-423CS**
- **Lockable cover - Item #F-720RLC**
- **Spine board assembly - Item #F-734RSA**
- **Gurney assembly - Item #F-730RSA**
- **Sling seat option - Item #F-706RLSS**
- **Wheelchair attachment - Item #F-705**
- **Transport cart - Item #F-714RTC**
- **Cycle attachment - Item #F-019CA**



Sling Seat Option



Wheelchair Option

The Revolution™ Lift is verified by an independent third party to be ADA compliant. This lift provides the ability to service both above-ground and in-ground pools and spas.

This extremely versatile lift is capable of providing access in a wide variety of situations.

The many adaptable components of The Revolution make it ideal for hospital and rehabilitation facilities.



U.S. Patent No.
D507,769 S

Item# F-700RL

FEATURES

- 500 lb (226.5 kg) weight capacity
- Independently verified to ADA Standards
- Rechargeable battery with wall mount charger
- Clears 30" high walls, up to 10" thick
- 360° powered rotation in either direction
- Adjustable footrest & lap belt
- Limited lifetime structural warranty
- 5 year pro-rated electronics warranty

PROUDLY SOLD BY:

The Revolution™ Lift

ADA Compliant