



Flame Red Frame  
Flame Blue Frame




**EASY HEIGHT ADJUSTMENT WITH SELF THREADING KNOB!**



10289RD



**COMES WITH CARRY BASKET AND TRAY!**



**BASKET CONVENIENTLY FOLDS FOR EASY STORAGE!**

## Deluxe 3 Wheel Aluminum Rollator, 8" Casters with Loop Locks

Available **1 Pack**

10289RD Flame Red, (Has replaced 10218), 1/cs  
10289BL Flame Blue, (Has replaced 10218), 1/cs  
815B Basket, 1/ea

### FEATURES:

- Comes standard with basket, tray and pouch.
- Ships in ultra compact carton. Comes with frame folded and compact. Simply loosen lock knob and pull trombone style frame to lock in place.
- Lightweight aluminum frame.
- Lightweight, solid 8" wheels for indoor or outdoor use.
- Adjustable handle height.
- **Special loop lock made of internal aluminum casting operates easily and ensures safety.**
- Easy, one hand folding.
- Limited Lifetime Warranty.




**EASY HEIGHT ADJUSTMENT WITH SELF THREADING KNOB!**



171



**COMES WITH CARRY BASKET AND TRAY!**



**BASKET CONVENIENTLY FOLDS FOR EASY STORAGE!**


## Deluxe 3 Wheel Steel Rollator, 8" Casters with Loop Locks

Available **1 Pack**

171 Chrome, (Has replaced 161), 1/cs  
815B Basket, 1/ea

### FEATURES:

- Comes standard with basket, tray and pouch.
- Ships in ultra compact carton. Comes with frame folded and compact. Simply loosen lock knob and pull trombone style frame to lock in place.
- Chrome plated steel frame.
- Lightweight, solid 8" wheels for indoor or outdoor use.
- Adjustable handle height.
- **Special loop lock made of internal aluminum casting operates easily and ensures safety.**
- Easy, one hand folding.
- Limited Lifetime Warranty.



**PARTS DIAGRAMS ON PAGE 33**

### SPECIFICATIONS

Item #	Handle Height (highest)	Handle Height (lowest)	Unit Width	Unit Length	Seat (depth)	Seat (width)	Seat (height)	Product Weight	Carton Shipping Weight	Weight Limit
<b>10289 SERIES</b>	38"	31"	25"	24"	NA	NA	NA	12 lbs.	20.5 lbs.	300 lbs.
171	38"	31"	25"	24"	NA	NA	NA	16.5 lbs.	20.5 lbs.	300 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.