

EASY HEIGHT ADJUSTMENT WITH SELF THREADING KNOB!



FOLD-UP & REMOVABLE BACKREST!



R726BL



SERRATED BRAKES!

PADDED SEAT OPENS TO A ROOMY CONVENIENT STORAGE POUCH! ALSO COMES WITH INSIDE ZIPPER POUCH ON UNDERSIDE OF PADDED SEAT!



COMES WITH NEW SEAMLESS PADDED SEAT WITH CARRY POUCH!



Blue

Red

Silver

Green

WEIGHS ONLY 12 LBS.

OPTIONAL BASKET AVAILABLE MODEL # 9501S72606

Aluminum Rollator with Fold Up and Removable Back Support, Padded Seat, 6" Casters and Loop Locks

Available **1 Pack**

- R726BL Blue, (Has replaced R6BL and 726BL), 1/cs
- R726RD Red, (Has replaced R6RD and 726RD), 1/cs
- R726SL Silver, (Has replaced R6SL and 726SL), 1/cs
- R726GR Green, 1/cs
- R726BK Black, 1/cs

FEATURES:

- Comes with new seamless padded seat with zippered pouch under seat.
- Removable, hinged, padded backrest can be folded up or down.
- 6" Black non-marring casters with soft grip tires are ideal for indoor and outdoor use.
- Comes with standard carry pouch under the seat for added privacy and security of personal items.
- Easy to use deluxe loop locks.
- Serrated brakes.
- Ergonomic handles are easy to grip and are adjustable in height.
- Basket sold separately (item # 9501S72606).
- Limited Lifetime Warranty.

Designer Winnie, Mimi Lite Deluxe Aluminum Rollator, Padded Seat, 6" Casters and Loop Locks

Available **1 Pack**

- 543 Flame Red Frame, 1/cs
- 544 Flame Blue Frame, 1/cs
- 548 Tan Frame, Tan Plaid Bag, 1/cs

FEATURES:

- Special loop lock made of internal aluminum casting operates easily and ensures safety.
- 6" casters with soft-grip tires are ideal for indoor and outdoor use.
- Tool Free removable padded backrest for comfort.
- Comes with standard carry pouch under the seat for added privacy and security of personal items.
- Comfortable padded seat.
- Handles are adjustable in height.
- Optional basket available (item # 830).
- Easy, one hand folding.
- Limited Lifetime Warranty.

EASY HEIGHT ADJUSTMENT WITH SELF THREADING KNOB!



548

SEE PAGE 21 & 22 FOR REPLACEMENT SEATS.

OPTIONAL BASKET AVAILABLE MODEL # 830


SEE PAGE 24



Flame Red Frame

Flame Blue Frame

Tan Frame



PARTS DIAGRAMS ON PAGE 29

SPECIFICATIONS

| Item # | Handle Height (highest) | Handle Height (lowest) | Unit Width | Unit Length | Seat Depth | Seat Width | Seat Height | Product Weight | Carton Shipping Weight | Weight Limit |
|-------------------------|-------------------------|------------------------|------------|-------------|------------|------------|-------------|----------------|------------------------|--------------|
| R726 SERIES | 37" | 32" | 24" | 28" | 14" | 12" | 21" | 12 lbs. | 16 lbs. | 300 lbs. |
| MIMI LITE DELUXE | 38" | 32" | 23" | 24" | 12" | 12" | 21.5" | 12 lbs. | 16 lbs. | 300 lbs. |

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

To order call toll free: 877.224.0946

ROLLATORS

drive™