

STRETCHAIR S675/S999



ISO 9001-2008 Certified

Our newest addition to the Bariatric Winco line offers versatility and safe transport

Stretchchair by Winco offers safe patient transport, lateral transfer aid and power vertical lift. Use as a rehabilitation chair, orthopedic chair, cardiac chair, geriatric chair or for bariatric patients.



KEY FEATURES

- Power Layflat positioning for transfer
- Flip up foot rest
- Independent leg rest
- Adjustable height foot rest
- Power Vertical lift
- Rechargeable, battery operated
- S675 – 675 lb. weight cap.
- S999 – 1000 lb. weight cap.
Power Recline and Transfer
and Vertical Lift

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

STRETCHAIR S675/S999

- 1 Conveniently positioned push handle
- 2 Stand assist grips for easier maneuverability
- 3 STANDARD IV pole mounts to right or left side
- 4 6" central locking & directional-locking casters
- 5 30"-38" vertical lift – layflat position
- 6 Adjustable height flip-up foot rests & independently adjustable leg rests
- 7 Wide, fully padded height adjustable arms
- 8 High resiliency foam for long-lasting comfort
- 9 Rechargeable battery operated hand control pendant converts chair to layflat position and operates the vertical lift.



OPTIONS



1B-133 TB



ADDITIONAL IV Pole IV

SPECIFICATIONS

S675

S999

Seat Height	22" (56cm)	22" (56cm)
Seat Width	28" (71cm)	31" (78.7cm)
Seat Depth	20" (51cm)	20" (51cm)
Overall Height	53.5" (136cm)	53.5" (136cm)
Overall Width	37" (94cm)	40" (101.6cm)
Back Height Above Seat	31.5" (80cm)	31.5" (80cm)
Overall Length (Upright)	45" (114.3cm)	45" (114.3cm)
Overall Length (Flat)	72" (or 76" w/ leg rests extended)	72" (or 76" w/ leg rests extended)
Adj. Arm Range or Fixed Height	30" to 40"	30" to 40"
Vertical Lift Range - Powered	Seat - 22" to 30", Layflat - 30" to 38"	Seat - 22" to 30", Layflat - 30" to 38" (700 Lb. Limit)
Weight	303 lbs. (137.4kg)	312 lbs. (141.5kg)
Weight Capacity	675 lbs. (306.2kg)	1000 lbs. (453.6kg) Power Recline and Transfer 1000 lbs. (453.6 kg) Vertical Lift


Range of Seat Angle stays parallel w/ floor

Range of Back Angle (Powered) 80 deg to zero (flat)

*Weight shown as a minimum-maximum range & varies depending on chair options ordered


STANDARD COLORS

Printed colors shown may vary slightly from actual color of product.



COLORS

OPTIONAL COLOR UPGRADES



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



11110111

REV D

009083