

Invacare eForms

For additional information, visit www.invacare.com/eForms

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order Literature: 800-828-6282

Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Please consider the environment before printing this page with your order form.
For more information on features and benefits of electronic order forms, please visit www.invacare.com/eForms

Invacare® Xtra™ - Ultralightweight Active Performance Folder

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



TRANSACTION TYPE:

- Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME TYPE

- XTRA Folding Square Frame with Swingaway Front End .. \$1,536.00 K0005

FRAME MODIFICATIONS

- HDPG Heavy Duty Package \$873.00 K0108

Heavy Duty Includes

- Weight Limitation 300 lbs.
- Reinforced seat and back upholstery
- Reinforced frame
- Triple-construction steel crossbraces

Width and Depth Matrix								
Depth Width	14"	15"	16"	17"	18"	19"	20"	22"
14"	X	X	X	X	X	S	S	S
15"	X	X	X	X	X	S	S	S
16"	X	X	X	X	X	S	S	S
17"	X	X	X	X	X	S	S	S
18"	X	X	X*	X*	X*	X*	X*	S
19"	S	S	X*	X*	X*	X*	X*	S
20"	S	S	X*	X*	X*	X*	X*	X*
22"	S	S	S	S	S	S	S	S

X - Sizes available without HDPG
 * - Sizes available with HDPG
 S - Sizes available through Specials

SEAT WIDTH

- W14 14" Wide \$45.00
- W15 15" Wide \$45.00
- W16 16" Wide No Charge
- W17 17" Wide \$45.00
- W18 18" Wide No Charge
- W19 19" Wide \$85.00
- W20 20" Wide \$300.00 E2201
- W22 22" Wide (Requires HDPG) \$300.00 E2201
- Only available with HDPG. Call Customer Service for special pricing and lead time for these sizes without HDPG.**

Arm Heights		
DESKADJCV, DESKADJSS, FULLADJCV, FULLADJSS	10" - 14"	Adjustable in 1" Increments
Fixed-Height Dual Point	10"	
Swingaway Padded Arm	8" - 10"	Adjustable in 1" Increments Add 1" with upholstered arm pads.
TAD Low Range	6" - 10"	Range adjusts in 1/2" increments.
TAF High Range	8" - 14"	Adjustable in 1" Increments

measurements from top of seat rail to top of armrest

SEAT DEPTH

- D14 14" Seat Depth \$45.00
- D15 15" Seat Depth \$45.00
- D16 16" Seat Depth No Charge
- D17 17" Seat Depth \$45.00
- D18 18" Seat Depth \$45.00
- D19 19" Seat Depth \$85.00
- D20 20" Seat Depth \$300.00 E2203
- D22 22" Seat Depth \$300.00 E2204
- Only available with HDPG. Call Customer Service for special pricing and lead time for these sizes without HDPG.**

ARM TYPE (PAIR)

- AT99 Swingaway Padded Arms No Charge
- DESKADJSS Spacesaver Adj Ht, Desk Length Arms \$195.00 E0973
Not available on 14" deep chairs.
- DESKFXDSS Spacesaver Fixed Ht, Desk Length Arms No Charge
Not available on 14" deep chairs. Not available with AAB1216 or AAB1720.
- FULLADJSS Spacesaver Adj Ht, Full Length Arms \$195.00 E0973
Not available on 14" deep chairs. Not available with AAB1216 or AAB1720.
- DESKADJCV Conv Adj Ht, Desk Length Arms \$195.00 E0973
Not available on 14" deep chairs.
- DESKFXDCV Conv Fixed Ht, Desk Length Arms No Charge
Not available on 14" deep chairs.
- FULLADJCV Conv Adj Ht, Full Length Arms \$195.00 E0973
Not available on 14" deep chairs.
- TAD T-Arm Adjustable Ht, Desk Length \$225.00 E0973
Bolt-On Arm Attachment.
- TAF T-Arm Adjustable Ht, Full Length \$225.00 E0973
Bolt-On Arm Attachment.
- LESSARMS Less Arms Less \$24.00

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® Xtra™ - Ultralightweight Active Performance Folder

T-ARM HEIGHT

Clear

- L6-10 Low Range 6" - 10" (T-Arm Height) No Charge
 H8-14 High Range 8" - 14" (T-Arm Height) No Charge

BACK TYPE

- A1418 10° Bend Back with Push Handles (14"-18") No Charge
 I014 Sportster Back without Push Handles (10"-14") \$30.00
 A1216 Straight Back with Push Handles (12"-16") \$30.00
 S1418 Straight Back with Push Handles (14"-18") \$30.00
 AAB1216 Adjustable Angle Back (12"-16") \$165.00 K0108
Not available with HDPG.
 AAB1720 Adjustable Angle Back (17"-20") \$165.00 K0108
Not available with HDPG.

BACK UPHOLSTERY TYPE

- SLIPBACK Slip On Back Upholstery No Charge
 ADJTENS Adjustable Tension Back Upholstery \$125.00 K0108
 L8601 Omit Back Upholstery Less \$10.00

ADJUSTABLE BACK TYPE SETTING

Clear

- ABH10 10" High Adjustable Back No Charge
 ABH11 11" High Adjustable Back No Charge
 ABH12 12" High Adjustable Back No Charge
 ABH13 13" High Adjustable Back No Charge
 ABH14 14" High Adjustable Back No Charge
 ABH15 15" High Adjustable Back No Charge
 ABH16 16" High Adjustable Back No Charge
 ABH17 17" High Adjustable Back No Charge
 ABH18 18" High Adjustable Back No Charge
 ABH19 19" High Adjustable Back No Charge
 ABH20 20" High Adjustable Back No Charge

CASTER OPTIONS (PAIR)

- 61DC 6" x 1" Urethane Casters No Charge
Not available with HDPG.
 1259 6" x 2" Semi. Pneumatic Casters \$60.00
 81DC 8" x 1" Urethane Casters No Charge
 815DC 8" x 1.5" Urethane Casters No Charge

Front Seat-to-Floor Heights				
	61DC	1259	81DC 815DC	815DC
6"	17.00			
	17 3/4			
	18 1/4			
	18 3/4			
	19.00			
8"	20.00			
		18.00	18 3/4	19.00
		19.00	19 3/4	20.00
		19 1/4	20.00	20 1/4
		20 1/4	21.00	21 1/4

Rear Seat-to-Floor Heights	
Wheel Size	Seat-to-Floor
20"	14 3/4" - 20"
22"	15 3/4" - 21"
24"	16 1/2" - 21 1/4"
Range adjusts in 7/8" increments.	

WHEEL OPTIONS (PAIR)

- 24COM7 24" Composite 7-Spoke Wheels Standard
 20COM 20" Composite Wheels No Charge
Not available with KIK tires.
 22COM 22" Composite Wheels No Charge
Not available with KIK tires.
 22SP 22" Spoke Wheels No Charge
Not available with HDPG.
Not available with KIK tires.
 24SP 24" Spoke Wheels No Charge
Not available with HDPG.
 L8701 Omit Rear Wheels Less \$50.00

TIRE OPTIONS (PAIR)

Clear

- PNEU Pneumatic Tires \$60.00
 PNEUFF Flat Free Inserts \$65.00 E2213
 URETH Urethane (Treaded) No Charge
Only available with composite wheels.
 KIKB KIK® Black Tires \$60.00

HANDRIM OPTIONS (PAIR)

- ALUMPC Aluminum Anodized Handrims Standard
 A1512 Black Plastic Coated Handrims \$72.00
 A1512K Black Plastic Coated Handrims \$72.00
Only available with 22SP and 24SP.
 AT1225-8 Plastic Coated with 8 Vertical Proj \$240.00
Only available with 22COM, 24COM7, and 24SP wheels.
 AT1224-8 Plastic Coated with 8 Oblique Proj No Charge
Only available with 24COM7 or 24SP.
 AT1225-12 Plastic Coated with 12 Vertical Proj \$240.00
Only available with 24COM7.
 L8702 Omit Handrims Less \$15.00

CAMBER ANGLE

Clear

- C00 0° Camber No Charge
 C03 3° Camber No Charge

MOUNTING HANDRIMS

Clear

- IN Mount Inside Handrim No Charge
Decrease overall width by 3/4"
 OUT Mount Outside Handrim Standard

WHEEL LOCKS (PAIR)

- A1598 Hill Holder Clamp-On Wheel Locks \$170.00 E0974
Not available with urethane tires or KIK black or pull-to-lock wheel locks.
 A1567 Pull to Lock Clamp-On Wheel Locks No Charge
 9691 Push to Lock Wheel Locks No Charge
Not available with AT5044 and ATL4A.

AXLE PLATE (PAIR)

Clear

- STDPT Standard Axle Plates Standard
 EXTPT Extended/Amputee Axle Plates \$72.00 E0959

AXLES

Quick release axles for spoked wheels are with dome button axle. Not available with HDPG.

- A1577 Quick Release Axles \$60.00
 PAXLE Permanent Axles No Charge
 A1573 Quad Release Axles \$89.00

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® Xtra™ - Ultralightweight Active Performance Folder

UPHOLSTERY AND COLORS

UPHOLSTERY COLORS

Monogramming and emblems available. Please check with Customer Service

- U240 Black Nylon Upholstery No Charge

FRAME COLORS

- 24P Wet Black No Charge
- 29P White Snow No Charge
- 60P Silver Metallic No Charge
- 115P Black Prism No Charge
- 105P Midnight Blue No Charge
- 116P Sunkissed Orange No Charge
- 117P Tangerine Red No Charge
- 30P Sunny Yellow No Charge
- 61P Electric Red No Charge
- 119P Bubblegum Pink No Charge
- 120P Grape Madness No Charge
- 62P Electric Blue No Charge
- 127P Grasshopper Green No Charge
- Touch Up Paint \$20.00

Touch up paint is not available in all colors.

Touch Up Color

Discontinued frame colors are available through our specials department.
C-070109030544

POSITIONING STRAPS AND ACCESSORIES

SEAT POSITIONING STRAPS (EACH)

- Clear
- A1515 Push Button Seat Positioning Strap \$50.00 E0978
 - A1321 Hook & Loop Seat Positioning Strap \$40.00 E0978

BACK ACCESSORIES

- BP02 Back Pack \$50.00
- 1483 Non-Adjustable Rigidizing Handle Bar \$85.00

SEAT ACCESSORIES

- L8651 Omit Seat Upholstery Less \$10.00
- SP-01 Seat Pouch \$35.00

ARM ACCESSORIES (PAIR)

Only available with AT99 or OMIT ARMS.

- A1571 Fabric Clothing Protectors \$46.00 K0108
- 1671-SM Rigid Side Guard, Small \$122.00 K0108
- 1671-LG Rigid Side Guard, Large \$122.00 K0108
- 1671-SMH Quad Rigid Side Guard, Small \$122.00 K0108
- 1671-LGH Quad Rigid Side Guard, Large \$122.00 K0108

MISCELLANEOUS ACCESSORIES

- 3421 3/4" to 1" Clamp Assembly (each) \$21.00
- A1209 Tire Pump..... \$32.00
- A1307 Spoke Protectors (for 20", 22", 24", 26")(pair) \$90.00 K0065
- Not available with projection handrims.
- A1569 Crutch & Cane Carrier \$46.00
- A1576 Wheel Lock Extensions 6" \$65.00 E0961
- A1580 Rear Anti-Tippers (pair) \$150.00 E0971
- A1575 Quick Release Caster Adapters (pair) \$125.00
- Not available with 8" casters.
- A1579 Front Anti-Tippers (pair) \$150.00 E0971
- Only available on tubular footplates.

FRONT RIGGINGS (PAIR)

Clear

- SWING1319 60° Swingaway Footrests No Charge
- SWING914 60° Swingaway Footrests MFX \$50.00
- Not available with longer pivot and slide tube option.
- 70NTAPER 70° Swingaway Footrests No Charge
- 70MFX 70° Swingaway Footrests MFX \$50.00
- Not available with longer pivot and slide tube option.
- 90SWING 90° Swingaway Footrests \$50.00
- Not available with longer pivot and slide tube option.
- AT5044 Elevating Legrests \$255.00 E0990
- AT5044 available in wet black only, must also use pull to lock wheel locks.
- ATL4A Center Pivot Articulating Legrests \$320.00 K0053
- Not available with longer pivot and slide tube options.
- 70STAPER 70° Tapered Footrests \$50.00
- 16"-20" wide frames only.

FOOTPLATE TYPE (PAIR)

Clear

- 1651 Flip-Up Composite Footplates No Charge
- Not available with 70MFX, 90SWING, and ATL4A.
- AT5543 Adjustable Angle Flip-Up Footplates \$165.00 K0040
- Dorsi/plantar angle adjustable.
- TUBULAR Individual Tubular Footplates \$59.00
- Not available with 70MFX, 90SWING, and ATL4A.
- A1564 Extra Large Aluminum Footplates \$59.00 K0041
- Only available with SWING1319, SWING914, AT5044, and ATL4A.
- A1562 Aluminum Footplates \$59.00
- Not available with with 70MFX, or 90SWING.

FOOTREST/LEGREST ACCESSORIES

- A1337 Calf Strap \$24.00 K0038
- IMPACTG Impact Guards \$20.00
- Standard on 70STAPER footrest. Not available with 90 degree, 70 degree MFX or 60 degree MFX swingaway footrests or elevating legrests.
- ALPT Longer Pivot and Slide Tubes \$84.00
- Composite or adjustable angle footplates only.
- A1600BK Heel Loop with Ankle Straps \$51.00 E0951
- A9342 "H" Strap \$51.00 K0039
- HLOOP Heel Loops without Ankle Straps \$51.00 E0951
- FRMPRB Chair Armor Blue Neoprene Frame Protector \$65.00
- FRMPRR Chair Armor Blue Neoprene Frame Protector \$65.00
- FRMPRS Chair Armor Blue Neoprene Frame Protector \$65.00
- FRMPRBK Chair Armor Black Neoprene Frame Protector ... \$65.00

Knee-to-Knee Range		
	Minimum	Maximum
70° Tapered	15 3/4"	19 3/4"
70°	13 1/2"	18 1/2"
70° MFX	8 1/2"	14 3/4"
Elevating	15 1/4"	21 1/4"
90°	5 1/4"	10 1/4"
60°	14 3/4"	18 1/2"
60° MFX	11"	15"
Range adjusts in 1" increments.		

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Xtra™ - Ultralightweight Active Performance Folder

SEATING CUSHION SELECTION Clear

For additional information (size availability, additional seating options, and pricing) see form #07-027 Rehab Seating and Positioning Products price list.

- EC Absolute Cushion \$104.00 E2601
Available in even sizes only. See form 07-027 for additional information.
- ITSG InTouch™ Stabilite™ Gentle Contour Cushion \$340.00 E2607
- ITSM InTouch Stabilite Max Contour Cushion \$340.00 E2605
- ITFG InTouch Flovair™ Gentle Contour Cushion \$360.00 E2607
- ITFM InTouch Flovair Max Contour Cushion \$360.00 E2607

SEATING BACK SELECTION Clear

For additional information (size availability, additional seating options, and pricing) see form #07-027 Rehab Seating and Positioning Products price list.

- PBR Personal Back 10™ Reg. Height \$565.00 E2615
- PBT Personal Back 10 Tall Height \$565.00 E2615
- PPR Personal Back 10 Plus Reg. Height \$645.00 E2615
- PPT Personal Back 10 Plus Tall Height \$645.00 E2615
- ITPR InTouch Propel™ Back \$438.00 E2611

Specifications	
Frame Style	Folding
Seat Width	14" - 20"
Seat Depth	14" - 18"
Overall Width Closed	12"
Adjustable Height Back	14" - 18"
Camber Angles	0°, 3°, 5°, 8°
Tubing O.D.	1" Aluminum
Weight Limit	250 lbs.
Heavy-Duty Package	300 lbs.
Weight of Chair	22 1/2 lbs
Warranty	Lifetime on frame and crossbraces

For additional seating needs please refer to the following:
Invacare Rehab Seating and Positioning Products Price List Form
#07-027.

Minimum Overall Width
Seat width plus 8 1/2" with 0° camber.
Seat width plus 10" with 3° camber.
*Measurements taken from outside of Handrim to outside of Handrim.

Weight Limitation 250 lbs.
Invacare recommends all Invacare wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

- Basic Chair Includes:**
- Folding Square Frame w/Swingaway Front Frame
 - Seat Widths (16" or 18" Wide)
 - Seat Depth (18" Deep)
 - Arm Type (AT99, DESKAFXDSS, DESKFXDCV)
 - A1418 10 Degree Bendback w/Push Handles (14" - 18")
 - Frame Finish
 - Slip on Back Upholstery
 - Caster Selection (61DC, 81DC, 815DC)
 - Fork Type
 - Wheel Options
 - URETH Urethane (treaded) Tires
 - ALUMPC Aluminum Anodized Handrims
 - Standard Axle Plates
 - Permanent Axles
 - Camber Selection (0 Degree or 3 Degree)
 - Front Rigging (SWING1319 or 70NTAPER)
 - Flip Up Composite Footplates
 - Knee to Heel Length Selection
 - ALPT Longer Pivot and Slide Tubes
 - Seat Upholstery
 - Wheel Locks (Push or Pull to Lock)

Comments

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.



Invacare Corporation
USA
One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

www.invacare.com

©2010 Invacare Corporation
All rights reserved. Trademarks are identified by the symbols ™ and ®.
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. KIK is a registered trademark of KIK Technology, Inc.
Specifications and pricing subject to change without notification.