

AQUATEC® *The* OCEAN

Aquatec Ocean joins the family of high-quality **Aquatec** products. The **Aquatec Ocean** is a height adjustable shower commode chair with 2 wheel options. **Ocean** comes standard with 5" locking swivel casters or is available as a self propel model with 24" wheels and push rims. The legs are adjustable, even above the 24" wheels. This provides more clearance over a toilet or for easier transfers. Height adjustment is easily done with push clips in 1" increments for a total of 5" of extra height.

Components

The **Ocean** is constructed from stainless steel for maintenance-free, rust-free use. The seat frame has a static 10° back tilt for comfort. The back upholstery is adjustable. Optional seat overlays are soft urethane. The seat, back upholstery and footrests are removable for cleaning. The entire chair can be disassembled for cleaning in a tub or clinical autoclave.

Features

- Stainless-steel frame
- Adjustable swing-away footrests
- Heel loops standard
- 5" casters with wheel locks standard
- 24" wheels have standard wheelchair brakes
- Can roll over most standard toilets
- Height adjustable flip-up arms
- Portable, disassembles for travel
- XL back provides additional 3" width and greater capacity

Warranty

Aquatec Ocean shower/commode chair frame is warranted for 2 years against manufacturer's defects, excluding upholstery and soft seats.



OceanXL arm detail



CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com

AQUATEC® *The OCEAN Accessories*



TOILET OVERLAY



SOLID INSERT



SPLASH GUARD



GEL SEAT



PADDED SEAT AND CHEST BELT



SAFETY SUPPORT BAR



SAFETY HARNESS COMBINATION



ARM LOCK



SIDE LATERAL AND ARM TROUGH



CALF OR AMP SUPPORT



CALF STRAP



CARRY BAG

SPECIFICATIONS	OCEAN (5" Casters)	OCEAN (24" Self Propel)	OCEAN XL
Width of chair	22.5"	27"	24.5"
Width between arms	18"	18"	21"
Height of arm, from seat	8.5" / 10"	8.5" / 10"	8.5" / 10"
Length of chair	37"	42.5"	37"
Height of chair	38" - 42.5"	38" - 42.5"	38" - 42.5"
Seat area	19.5" x 18"	19.5" x 18"	19.5" x 18"
Seat height	19" - 23"	19" - 23"	19" x 23"
Standard wheel diameter	5"	5"	5"
Accessory wheel diameter	24"	24"	24"
Height adjustment increments	1"	1"	1"
Distance between floor and seat frame	16" min. to 21" at max. leg height	16" min. to 21" at max. leg height	16" min. to 21" at max. leg height
Min. width between side sections, under frame, for toilet clearance	15"	15"	15"
Capacity	290 lb.	290 lb.	400 lb.
Weight of chair	33.5 lb.	33.5 lb.	33.5 lb.
Back support	16"W x 14"H	16"W x 14"H	16"W x 14"H
Standard hygiene seat opening	8" wide x 15.5" long	8" wide x 15.5" long	8" wide x 15.5" long
Product number	A710900	A712900	A713900

CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com