



ECLIPSE 2010 ORDER FORM

COLOURS



tel: 951-808-9131 Fax: 951-808-9949
Toll Free: 800-892-8998 www.colourswheelchair.com

BILL TO

Company	P.O.#	Date	S/N	for Office Use Only
Completed by	Address			
Email Address				
City	State	Zip	Phone	Fax
Disability	Weight	Height		
Customer Name	Client Signature X	Email		

SHIP TO

Name	Company	Phone
Address	City	State Zip

A MEASUREMENTS

FRAME No camber frame

- ☐ (standard no-camber frame with tapered footrest) \$ std
- ☐ (custom inside foot width _____) or ☐ no taper frame \$ 267
- ☐ (axle plates mounted inside -reduces overall width by 3.0") \$ 112
- Seat Width 13 14 15 16 17 18 \$ std
1. Seat Width **E2201** 19 20 \$ 300
- Seat Depth 12 13 14 15 16 (Short Frame) \$ std
2. Seat Depth 14 15 16 17 18 19(3" Longer Frame) \$ std
- Seat Depth **E2203** 20 (Custom Frame) \$ 250
- BACKREST** (adj. range) (adj. range)
- Height 8 | 9 10 11 12 | 13 14 15 | \$ std
3. 16 | 17 | 18
- Type ☐ Folding with adjust. height and angle \$ std
- ☐ NEW Locking mechan. (adjust. angle) \$ std

AXLE PLATE

- ☐ One Slot, Reversible \$ std
- ☐ One Slot, Amputee \$ 107
- ☐ Dual Receivers \$ 107

Camber Front Receiver (0°-14°) _____°

Camber Rear Receiver (0°-14°) _____°

FOOT PLATES

- ☐ Standard single loop (with riser if necessary) \$ std
- ☐ Single Flip-Back - single loop (no riser avail.) **K0040** \$ 150

Base Wheelchair Price (K0005) \$ 2,241

Total Section A \$ _____

Total Section B \$ _____

Total Section C \$ _____


Total Section D \$ _____

Custom Modifications Quoted \$ _____

Order Grand Total \$ _____

- A. Cushion Thickness _____" 1" - 4"
- B. Camber Angle _____° (0°-14°)
- C. Lower Leg Length _____"
- D. Inside Wheel Clearance _____" 1" is std.
4. Front Seat Height _____" 17.5" - 19.5" (w/3" caster)
5. Rear Seat Height _____" 13.5" - 16.0" (w/24" rear wheel)
6. Overall Width _____"
- Std. is 2.5"
8. Center of Gravity _____" (subtract 5" off seat width w/tapered)
9. Footrest Width _____" frame to get inside foot clearance)

Custom Dimensions Available. Contact Colours For Details.



NOTES

B WHEELS & TIRES

- ☐ Standard Radial spoke wheel ☐ 24" ☐ Grey Primo hi-press. tire ☐ Standard low press. tire \$ std
- ☐ Heat-Treated wheel * ☐ 24" ☐ 25" ☐ 26" \$ 252
- ☐ Cross Knot Laced Spokes (sport)* ☐ 24" ☐ 25" ☐ 26" ☐ 27" \$ 211
- ☐ Twister Super Strength wheel * ☐ 24" ☐ 25" ☐ 26" \$ 310
- ☐ Spinergy LX 12 spokes everyday * ☐ 24" ☐ 25" ☐ 26" Spokes Colors ☐ Black ☐ Yellow ☐ Blue ☐ Red \$ 895
- ☐ Flex Handrim (available w/Spinergy 24" & 25" 12 spoke only) ☐ 24" or ☐ 25" add \$ 450
- ☐ Bite Laced Wheels with NEW locking hubs (patent pending) ☐ 24" or ☐ 25" with h.p. tire \$ 895
- ☐ X-Core (3-spoke comp.)* ☐ 25" ☐ Black(std.) ☐ Yellow ☐ Blue ☐ Red ☐ White (add \$200) camber 0-2° max \$ 743
- ☐ Y-360 Billet Wheels (available in 25" only) ☐ Saw Blade ☐ Hearts ☐ T-Style camber 0-2° max \$ 999
- ☐ MOD (\$1,499) ☐ Tribal Sun ☐ Girls ☐ Eagle ☐ Holy Cross (Specify color) Y-360 Spinners ☐ Clover (Add) \$ 650
- ☐ Y-360 Billet Wheels with Bite Locking Hub (available in for all styles) ☐ Mod-Style add \$600 \$ 1600
- ☐ E2213 Flat-Free Inserts (24" only) ☐ 24" \$ 117
- ☐ Solid tires (dark grey) ☐ 24" ☐ 25" ☐ 26" \$ 135
- ☐ E2211 Marathon Schwalbe tires (High Pressure) ☐ 24" ☐ 25" (All Black Tires) \$ 99
- ☐ E2211 Knobby ☐ 24x1.375 mod. tread ☐ 24x2.0 aggr. tread ☐ 25x1.95 aggr. tread \$ 117

* Standard with black wall hi-press. tire and anodized handrims

Refer to Section D for more options!

Section B Total \$ _____

Valid From Aug '10 . Prices &
Specifications Are Subject To Change
Without Notice.

ECLIPSE '2010

COMPONENTS & ACCESSORIES

C

Handrims

<input type="checkbox"/> Standard For 24", 25", or 26" Wheel	Anodized Black	\$ std
<input type="checkbox"/> Vinyl Coated	Black	\$ 142
<input type="checkbox"/> The Natural-Fit Handrim <input type="checkbox"/> 24" or <input type="checkbox"/> 25" with <input type="checkbox"/> reg. grip or <input type="checkbox"/> super grip +\$55		\$ 295
<input type="checkbox"/> The Natural-Fit Handrim (LT) <input type="checkbox"/> 24" or <input type="checkbox"/> 25" with <input type="checkbox"/> reg. grip or <input type="checkbox"/> super grip +\$85		\$ 325

Forks

<input type="checkbox"/> E1015 Frog Leg Suspension System (Std. 1.0 Caster)	<input type="checkbox"/> 3" <input type="checkbox"/> 4" <input type="checkbox"/> 5" <input type="checkbox"/> Black <input type="checkbox"/> Silver	\$ 535
<input type="checkbox"/> E1015 Frog Leg Suspension System (Wide Body 1.5 Caster)	<input type="checkbox"/> 4" <input type="checkbox"/> 5" <input type="checkbox"/> Black <input type="checkbox"/> Silver	\$ 535
<input type="checkbox"/> E1015 Suspension Springs, Front Casters (may increase front seat height 1")		\$ 55

Caster Options (installed on new chair ordered)

Standard Casters <input type="checkbox"/> 3" Micro Caster (black, clear, red) _____ <input type="checkbox"/> 4" Solid Caster <input type="checkbox"/> 5" Solid Caster	\$ std
Lighted Casters <input type="checkbox"/> 3" Lighted <input type="checkbox"/> 4" Lighted <input type="checkbox"/> 5" Lighted	\$ 65
<input type="checkbox"/> 1.0 Wide Caster Wheels Alum. Hub (pr) <input type="checkbox"/> 3" <input type="checkbox"/> 4" <input type="checkbox"/> 5" <input type="checkbox"/> Black Hub <input type="checkbox"/> Silver Hub	\$ 160
<input type="checkbox"/> 1.5 Wide Caster Wheels Alum. Hub (pr) <input type="checkbox"/> 4" <input type="checkbox"/> 5" <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> Yellow	\$ 160
<input type="checkbox"/> 1.5 Wide Caster Wheels Composite Hub (pr) <input type="checkbox"/> 4" <input type="checkbox"/> 5" <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> Yellow	\$ 120

Push Handles

<input type="checkbox"/> Integrated w/ Fixed Backrest Height	\$ 126
<input type="checkbox"/> Bolt-On <input type="checkbox"/> Quick Release (Add \$48) specify ht. from floor _____"	\$ 142
<input type="checkbox"/> Slider Type Quick Release Overall Length <input type="checkbox"/> 12" <input type="checkbox"/> 16" <input type="checkbox"/> 20"	\$ 142
<input type="checkbox"/> Adjustable Height Push Handles (not available with armrest)	\$ 199

Arm Rests

<input type="checkbox"/> E0973 Aluminum, Padded - Swing Away Height Adjustable	\$ 225
--	--------

Miscellaneous

<input type="checkbox"/> E0971 Anti-Tippers (pair) Quick Release	\$ 112
<input type="checkbox"/> PRO-TEX Fender Guards - Black ABS (Bolt- on only) <input type="checkbox"/> Cut to fit Specs. <input type="checkbox"/> Do not Install	\$ 175
<input type="checkbox"/> Rigid Side Guards - Black ABS <input type="checkbox"/> removable <input type="checkbox"/> non-removable <input type="checkbox"/> cut to wheel (add \$98)	\$ 142
<input type="checkbox"/> Rigid Side Guards - Alum. <input type="checkbox"/> removable <input type="checkbox"/> bolt-on <input type="checkbox"/> pol. alum. <input type="checkbox"/> Black	\$ 197
<input type="checkbox"/> Fabric Clothing Guards Black Nylon w/Velcro	\$ 78
<input type="checkbox"/> Frame Protectors Black Leather/Zipper type	\$ 78
<input type="checkbox"/> K0065 Spoke Protectors <input type="checkbox"/> 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26" Black	\$ 107
<input type="checkbox"/> E0978 Seat Belt w/Auto Buckle 1 1/2" Belt	\$ 69
<input type="checkbox"/> K0038 Calf Strap - adjustable length Black	\$ 24

Frame Colors (Powder-coated)

<input type="checkbox"/> Gloss Black	<input type="checkbox"/> Candy Blue	<input type="checkbox"/> White	<input type="checkbox"/> Yellow	<input type="checkbox"/> Cobalt Blue	<input type="checkbox"/> Candy Red
<input type="checkbox"/> Twilight Flake	<input type="checkbox"/> Candy Purple	<input type="checkbox"/> Candy Green	<input type="checkbox"/> Vermilion	<input type="checkbox"/> Fire Engine Red	<input type="checkbox"/> Reed Green
<input type="checkbox"/> Neon Red	<input type="checkbox"/> Neon Green	<input type="checkbox"/> Orange Sparkle	<input type="checkbox"/> Magenta	<input type="checkbox"/> Candy Raspberry	<input type="checkbox"/> Mellon Yellow
<input type="checkbox"/> Pigeon Blue	<input type="checkbox"/> Matte Grey	<input type="checkbox"/> Matte Purple	<input type="checkbox"/> Matte Blue	<input type="checkbox"/> Gun Smoke Silver	<input type="checkbox"/> Blue Lilac
<input type="checkbox"/> Matte Green	<input type="checkbox"/> Matte Red	<input type="checkbox"/> Matte Black	<input type="checkbox"/> Copper	<input type="checkbox"/> Sparkle Sky Blue	\$ std
<input type="checkbox"/> Granny Sparkle Green					

Custom Color (Custom Polished)

<input type="checkbox"/> Polished Alum. Frame Only	\$ 800
--	--------

Upholstery Choice & Color

<input type="checkbox"/> Standard Adjustable Ht. Black	<input type="checkbox"/> Black	\$ std
<input type="checkbox"/> K0108 Adjustable Tension Backrest	<input type="checkbox"/> Black	\$ 135
<input type="checkbox"/> Custom Colored (for std. or adjust. style)	<input type="checkbox"/> Red <input type="checkbox"/> Blue	\$ 60
<input type="checkbox"/> Fabrics for FUN (for std. or adjust. style)	<input type="checkbox"/> 001 <input type="checkbox"/> 002 <input type="checkbox"/> 003 <input type="checkbox"/> 004 <input type="checkbox"/> 006 <input type="checkbox"/> 010 <input type="checkbox"/> 011	
	<input type="checkbox"/> 014 <input type="checkbox"/> 015 <input type="checkbox"/> 016 <input type="checkbox"/> 017 <input type="checkbox"/> 018 <input type="checkbox"/> 020 <input type="checkbox"/> 021	
	<input type="checkbox"/> 022 <input type="checkbox"/> 023 <input type="checkbox"/> 024 <input type="checkbox"/> 025 <input type="checkbox"/> 026 <input type="checkbox"/> 027 <input type="checkbox"/> 028	\$ 95
<input type="checkbox"/> Backrest Leather Trim (outside panels)		\$ 77

* <input type="checkbox"/> RazorBack (Carbon Backrest) By ADI	<input type="checkbox"/> TL-Style 16" fixed (13-18 wide) <input type="checkbox"/> Adj. Model (add \$150)	\$ 700
	<input type="checkbox"/> Ls-Style 10" fixed (13-19 wide) <input type="checkbox"/> Adj. Model (add \$150)	\$ 600

* No secondary discount on Price

☛ May not fit and function on all chairs depending on wheelchair configuration.

Section C Total

\$ _____

You're only here for a short visit. Don't hurry and don't worry. And make sure you smell the flowers along the way.





Colours in Motion warrants all rigid hardware for the lifetime of its registered owner. (Please see warranty registration card for details.)



D

EXTRA OPTIONS FOR ALL MODELS

Seat Cushion

<input type="checkbox"/> Foam Cushion w/ Nylon Cover	<input type="checkbox"/> 1" <input type="checkbox"/> 2" <input type="checkbox"/> 3" <input type="checkbox"/> 4"	\$ 78
<input type="checkbox"/> The Curve Cushion (<i>Rehab Co.</i>)	<input type="checkbox"/> 3" 	\$ 133
<input type="checkbox"/> E0734 ROHO Cushion Quadtro	<input type="checkbox"/> 2" Low Profile  <input type="checkbox"/> 4" High Profile  (USA only)	\$ 525
<input type="checkbox"/> E0734 ROHO Nexus Spirit Cushion	<input type="checkbox"/> 2" Low Profile  (USA only)	\$ 525

Other (installed on new chair ordered)

<input type="checkbox"/> Bicycle Plates with extra pair of axle receivers for Spazz, Blade, Chump and Zephyr	\$ 250
<input type="checkbox"/> Wheelchair Tray (size avail.) <input type="checkbox"/> 11" x 13" <input type="checkbox"/> 13" x 15" <input type="checkbox"/> 13" x 17" <input type="checkbox"/> Raw <input type="checkbox"/> Black <input type="checkbox"/> Polish Alum. (Add \$250)	\$ 250
<input type="checkbox"/> Ski Tips on Caster Forks	\$ 275
<input type="checkbox"/> AZ-1 Portable Hand Controls	\$ 349
<input type="checkbox"/> Anti Skid bar for rear of <input type="checkbox"/> Boing! <input type="checkbox"/> Spazz <input type="checkbox"/> Krypto (not available w/anti-tippers)	\$ 99
<input type="checkbox"/> Basketball Roller	\$ 62
<input type="checkbox"/> One Arm Drive Custom Wheelchair (<i>Zephyr, Spazz, Blade, Hammer, Chump</i>)	\$ CFP
<input type="checkbox"/> Heavy Duty Wheelchair for 300 or 350 lbs. (<i>Spazz, Avenger, Eclipse</i>)	\$ CFP
<input type="checkbox"/> Tie Down Brackets 4pcs. <i>Spazz, Chump, Spazz-G, Chump-G (Have been crash tested wc-19)</i>	\$ 210

Custom Paint Scheme (installed on new chair ordered)

<input type="checkbox"/> Polish All Parts excluding the frame	\$ CFP
<input type="checkbox"/> Powder Coat bolt on parts	\$ CFP
<input type="checkbox"/> 24 karat gold plating on various parts	\$ CFP

Matching Gear (provide extras - installed on chair)

<input type="checkbox"/> Matching Fabric on Leg Strap - select fabric _____	\$ 55
<input type="checkbox"/> Matching Fabric on Padded Side Guards - select fabric _____	\$ 210
<input type="checkbox"/> Cushion Cover (<i>Std. Cushion Only</i>) - select fabric _____ Width _____ Depth _____ Thickness _____	\$ 90

Spare Parts (provide extras - not installed on chair)

<input type="checkbox"/> Colored Tires (<i>pr</i>) <input type="checkbox"/> 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26" <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> Yellow	\$ 65
<input type="checkbox"/> Marathon Schwalbe Tires (<i>pr</i>) <input type="checkbox"/> 24" <input type="checkbox"/> 25" <input type="checkbox"/> Black	\$ 80
<input type="checkbox"/> Knobby Tires/Tubes (<i>pr</i>) <input type="checkbox"/> 24x1.375 mod. tread (<i>Gray</i>)	
<input type="checkbox"/> 24x2.0 aggr. tread (<i>Gray with <input type="checkbox"/> Blue or <input type="checkbox"/> Red Wall</i>) <input type="checkbox"/> 25x1.95 aggr. tread (<i>All Black</i>)	\$ 120
<input type="checkbox"/> High Pressure Tubes (<i>pr</i>) <input type="checkbox"/> 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26" <input type="checkbox"/> 27"	\$ 12
<input type="checkbox"/> Bearings Rear Wheels or Caster Housing (<i>set of 4</i>)	\$ 12
<input type="checkbox"/> Bearings for Caster Wheels (<i>set of 4</i>)	\$ 10
<input type="checkbox"/> Quick Release Big Button Rear Axles (<i>pr.</i>) <input type="checkbox"/> 4.5" <input type="checkbox"/> 5.5"	\$ 60
<input type="checkbox"/> K0108 Quad Wheel Quick Release Flip Up	\$ 99
<input type="checkbox"/> Front Axles <input type="checkbox"/> Front - quick release	\$ 107

Merchandise (provide extras - not installed on chair)

<input type="checkbox"/> Women S/S Tee <input type="checkbox"/> S <input type="checkbox"/> M <input type="checkbox"/> L <input type="checkbox"/> XL <input type="checkbox"/> 2XL <input type="checkbox"/> Black <input type="checkbox"/> White <input type="checkbox"/> Light Blue	\$ 30
<input type="checkbox"/> Mens S/S Tee <input type="checkbox"/> S <input type="checkbox"/> M <input type="checkbox"/> L <input type="checkbox"/> XL <input type="checkbox"/> 2XL <input type="checkbox"/> Black <input type="checkbox"/> White <input type="checkbox"/> Gray	\$ 30
<input type="checkbox"/> Caps Flex Fit <input type="checkbox"/> S/M <input type="checkbox"/> L/XL <input type="checkbox"/> Black <input type="checkbox"/> Gray	\$ 24
<input type="checkbox"/> Beanies <input type="checkbox"/> one size fits all <input type="checkbox"/> Black <input type="checkbox"/> Gray	\$ 20
<input type="checkbox"/> Pillows with insert 14" x 14" <input type="checkbox"/> Blk Flame <input type="checkbox"/> Camo Girls <input type="checkbox"/> Corazones <input type="checkbox"/> High Speed	\$ 55

Gear (provide extras - not installed on chair)

<input type="checkbox"/> Grab Bar Protector Removable (<i>3/4" backrest tubing</i>)	\$ 35
<input type="checkbox"/> Wheel Bag <input type="checkbox"/> 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26" <input type="checkbox"/> 27"	\$ 350
<input type="checkbox"/> Colours Backpack (<i>one size only</i>) <input type="checkbox"/> Black w/red logo <input type="checkbox"/> Black w/blue logo	\$ 55
<input type="checkbox"/> Positioning strap <input type="checkbox"/> Velcro one end <input type="checkbox"/> Velcro both ends _____ Length	\$ 85
<input type="checkbox"/> Ratchet Strap <input type="checkbox"/> small <input type="checkbox"/> medium <input type="checkbox"/> large	\$ 125

Wheelchair Fingerless Gloves (*pr.*)

<input type="checkbox"/> Sport Grip <input type="checkbox"/> Exec. Style <input type="checkbox"/> Smooth Style	\$ 29
<input type="checkbox"/> S <input type="checkbox"/> M <input type="checkbox"/> L	

Section D Total \$ _____

CFP = Call For Price

 May not fit and function on all chairs depending on wheelchair configuration.

































You're only here for a short visit. Don't hurry and don't worry. And make sure you smell the flowers along the way.

Colours in Motion warrants all rigid hardware for the lifetime of its registered owner. (Please see warranty registration card for details.)

C O L O R S & F A B R I C S

A C O L O R S

Valid From Aug '10. Prices & Specifications Are Subject To Change Without Notice.

<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>	
	2 1 0 0 gloss black		2 1 0 2 candy blue		2 1 0 6 white		2 1 0 9 yellow		2 1 4 7 orange sparkle		2 1 1 1 candy red		2 1 5 1 cobalt blue
<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>	
	2 1 1 3 candy purple		2 1 3 0 candy green		2 1 4 8 noen red		2 1 0 5 fire engine red		2 1 4 9 noen green		2 1 5 0 vermillion		2 1 3 8 gun smoke silver
<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>	
	2 1 5 2 melon yellow		2 1 3 1 magenta		2 1 5 3 reed green		2 1 5 4 pigeon blue		2 1 5 6 blue lilac		2 1 5 7 sparkle green		2 1 5 8 candy raspberry
<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>	
	2 1 5 0 matte gray		2 1 5 1 matte purple		2 1 5 2 matte blue		2 1 5 3 matte green		2 1 5 4 matte red		2 1 5 5 matte black		2 1 5 9 Sparkle Sky Blue
<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>							
	2 1 1 2 twilight flakes		6 0 0 0 polish alum.		6 0 0 1 gold plate		2 1 5 9 copper						

B F A B R I C S

<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>	
	F - 0 0 1 Black Flames		F - 0 0 2 Camouflage Girls		F - 0 0 3 Claws		F - 0 0 4 Corazones		F - 0 1 0 High Speed		F - 0 1 1 Skulls On Fire		F - 0 1 4 Skate Line
<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>	
	F - 0 1 5 Flame Rider		F - 0 1 6 Girls On Fire		F - 0 1 8 Vegas		F - 0 1 7 Dragon Tattoo		F - 0 0 6 Wheels On Fire		F - 0 2 0 Wild Thing		F - 0 2 1 Bloom
<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>	
	F - 0 2 2 Curve Ball		F - 0 2 3 Pink Lily		F - 0 2 4 Pulling G's		F - 0 2 5 Ruff Rider		F - 0 2 6 Little Scouts		F - 0 2 7 Virgin Mary		F - 0 2 8 Lady In Red