


- Designed for doorways that swing in or out
- Adjusts to a range of threshold rises for a perfect fit
- Durable welded construction
- Slip resistant grooved aluminum
- Accommodates wheelchairs and scooters with various wheel configurations
- Aligns directly against the door threshold up to 6"



Range of Threshold Heights
All In Each Ramp!



The Ramp with Infinite Heights

LIMITED **LIFETIME**
WARRANTY

Item Number	Length	Width	Threshold Height	Weight	Capacity*
ATH1232	12"	32"	1" - 2"	6 lbs.	800 lbs.
ATH1236	12"	36"	1" - 2"	12 lbs.	800 lbs.
ATH2432	24"	32"	2" - 3½"	12 lbs.	800 lbs.
ATH2436	24"	36"	2" - 3½"	12 lbs.	800 lbs.
ATH3632	36"	32"	3½" - 5"	12 lbs.	800 lbs.
ATH3636	36"	36"	3½" - 5"	12 lbs.	800 lbs.
ATH4832	48"	32"	4" - 6"	12 lbs.	800 lbs.
ATH4836	48"	36"	4" - 6"	12 lbs.	800 lbs.

* Maximum ramp capacity: 400 lbs. (one axle) – 800 lbs. (two axles).