

S-305

Ergonomic Class

ULTRA LIGHTWEIGHT MANUAL WHEELCHAIR

SPECIFICATIONS	S-305-16	S-305-18
Seat Width	16 in.	18 in.
Seat Depth	17 in.	17 in.
Armrest Height	9-11 in.	9-11 in.
Seat Height	19 in.	19 in.
Back Height	17 in.	17 in.
Overall Height	38 in.	38 in.
Overall Open Width	24 in.	26 in.
Folded Width	12 in.	12 in.
Overall Length	43 in.	43 in.
Weight w/o Riggings	29 in.	29 in.
Weight Capacity	250 lbs.	250 lbs.
HCPCS Code	K0005	K0005

(SUBJECT TO CHANGE WITHOUT NOTICE)



ULTRA LIGHTWEIGHT **SPECIALIZED** FRAME...

- 24" Quick Release rear spoke/polyurethane/high profile/flat free wheels
- Axle Adjustable Seat Height (18", 19", 20")
- One piece foldable footplate
- Folding backrest, for transporting chair
- Swing in & away footrest for maximum safety while entering or exiting the chair
- S-Shape Frame (patented)
- Anti-Bacterial Upholstery/Cushion
- Frame Color: Pearl Silver & Rose Red
- Seat Width: 18"x17" OR 16" x 17"