

P1 - Wednesday, April 7, 2021 7:57:40 AM

PRODUCT NAME: Admiral[™] Pool/Spa Access Lift

CSI SECTION: 13 11 46 Swimming Pool Accessories **APPLICATIONS:**

In Ground Pools & Spas with flush pool decks having up to **2-inch (5.08cm)** high curbs. •

SUBMITTALS:

 Alternate make or model lifts will not be considered unless equal to the specified product in every respect. Alternative products must be submitted for approval no less than ten (10) calendar days prior to the bid due date. Submittal data must include complete documentation relating to all the specified features and include all manufacturer's sales literature, specification sheets, installation/operation manuals, engineered drawings and documents that support the strength, service capacity, and the ability to meet the appropriate Federal Regulatory Standards associated with the product.

QUALITY ASSURANCE:


- The product must be a standard product of a company registered with the United States Food and Drug Administration (FDA) as a class 1 medical device manufacturer and that is regularly engaged in the manufacture of commercial grade swimming pool and spa accessibility equipment.
- The product shall be structurally capable of supporting a **450 pound (204.12kg)** load capacity and static load tested to 1.5 times the rated load.
- The product shall be certified to be compliant with UL test standard 60335-2-1000 Standard of Household and Similar Appliances: Particularly Requirements for Electrically Powered Pool Lifts.
- The product shall be tested and confirmed to comply with design standards outlined in the United States ADA Accessibility Standards for Accessible Designs.

SPECIFICATIONS:

- The products main structure shall be constructed of type 304 stainless steel; no dissimilar metallic components shall be allowed. The product shall be finished with a minimum 4-mil thick, UV resistant powder coat finish. All fastening hardware shall be of type 316 stainless steel. Product shall be supplied with a Linak Actuator suitable for medical use, a separate wall mount battery charger, a hand-held controller, and a 24-volt rechargeable battery.
- The chair of the product shall be of a rigid one-piece, polyethylene, molded design. The minimum seat width shall be **18-inches (45.72cm)** and the seat back shall be a minimum height of **17.5-inches** (44.45cm). The seat shall feature dual flip-up arms and supplied with a 2-inch (5.08cm) wide adjustable safety belt with hook and loop fasteners.
- The seat height of the product will be adjustable from a load height of 18-1/4-inches (46.36cm) to 24-1/4inches (61.60cm) in 2-inch (5.08cm) increments to accommodate different height wheelchairs.
- The product shall be capable of accommodating a deck to water distance (water draft) of up to 8-inches (20.32cm) and capable of an anchor set back distance of up to 19-inches (48.26cm).
- Product can be field reversible and includes an adjustable and removable footrest that moves with the seat.

WARRANTY:

Product shall have a minimum 5 year structural and 5 year pro-rated electronic warranty.



LIFT PART NUMBERS:

ltem No	Order No	Description
F-ADMRL	INV01513	Admiral Lift White Powder
F-ADMRL-T	INV01757	Admiral Lift White Powder
F-ADMRL-G	INV01756	Admiral Lift White Powder
F-ADMRL-W	INV01755	Admiral Lift White Powder
F-ADMRL-C	INV01907	Admiral Lift (Specify Seat ar

ANCHORING OPTIONS:

ltem No	Order No	Description
F-04CAJP	INV01411	Admiral Lift Anchor for Star
F-08CAJP	INV01414	Admiral Lift Anchor for Pav
F-046QAB-AT1	INV01406	Admiral Lift Anchor for Star
F-40220-AT1	INV01444	Admiral Lift Anchor for Wo

PRODUCT OPTIONS:

Item No	Order No	Description
F-0904CVR	INV01990	Admiral Lift Protective Sea
F-120PPC	INV00525	Admiral Lift Protective Lift
F-120PPC-T	INV00528	Admiral Lift Protective Lift
F-120PPC-G	INV00527	Admiral Lift Protective Lift
F-044SCH	INV01693	Admiral Lift Solar Charger
F-422PLH-2	INV01455	Admiral Lift Adjustable Hea
F-423CS	INV01457	Admiral Lift Chest Strap
F-428HT	INV01458	Admiral Lift Transport Cart
F-LRB	INV01422	Admiral Lift Pull Out Leg Re
F-LRT	INV01952	Admiral Lift Pull Out Leg Re
F-LRG	INV01954	Admiral Lift Pull Out Leg Re
F-LRW	INV01953	Admiral Lift Pull Out Leg Re
F-019CA	INV01379	Admiral Lift Cycle Attachm

Aqua Creek Products LLC 9889 Garrymore Lane Missoula, MT 59808 www.aquacreekproducts.com Tel: 888-687-3552 / 406-549-0769

Fax: 406-541-8958



9889 GARRYMORE LANE MISSOULA, MT 59808 (406) 549-0769 FAX (406) 549-2602 www.AguaCreek.com

PROPRIETY AND CONFIDENTIAL INFORMATION

THIS DRAWING AND ALL INFORMATION CONTAINED HEREIN IS THE SOLE PROPERTY OF AQUA CREEK PRODUCTS LLC. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT WRITTEN PERMISSION FROM AQUA CREEK IS PROHIBITED.

***SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE**

Coat w/Blue Seat	
Coat w/Tan Seat	
Coat w/Gray Seat	
Coat w/White Seat	
nd Powder Coat Colors)	

indard Concrete Decks ver Applications andard Concrete Decks with Quick Attach Knobs ood Deck Applications

at Cover in Blue Cover in Blue Cover in Tan

Cover in Grav

ead Rest

lest Blue Rest Tan Rest Gray

Rest White

ient



