



The MAXI+

Maxi+ is designed for the heavier, larger or taller user. The wider frame provides more seat room and a larger walking area. Dolomite has the strongest frames on the market; **Maxi+** is rated for a 440 lb. capacity. **Maxi+** standard size is the tallest model in the **Dolomite** line. Choose **Maxi+** for a stable, walking experience and travel safely in style.

Components

Maxi+ is made from high tensile steel with welded joints for strength and low maintenance. **Maxi+** is comes in Metallic Blue epoxy-based powder coat finish. **Dolomite Maxi+** comes in two frame heights, see chart for specs. The frames fold front to back for storage. A safety latch prevents accidental folding when using the seat handle.

Maxi+ includes a wide flip-up seat for more walking room. Anatomical grips have a larger surface which reduces hand fatigue for arthritic users. Flexible back strap is adjustable. Large 8" wheels roll easily indoors or outdoors. Red Curb climber on the rear wheel assists when lifting the walker over a curb or step.

Braking

Maxi+ comes with the patented DABS braking system, with adjustable tension for varied grip strength. Brake handles have a wide flange and large loops for comfortable hand positioning. The brake mechanism is located within the fork for protection. When parked and using the seat, pushdown the handles to set the parking brakes.

Features

- Extra wide, 20.5" between handles
- Large flip-up seat
- 440 pound capacity
- Perfect for the taller person, 25.5" seat height
- Large 8" wheels
- Handles adjust 7" in 1" increments
- Large basket for personal items
- Adjustable back support

Warranty

Dolomite warrants the **Maxi+** frame to be free from manufacturer's defects for 2 years. The DABS brake handles have a 5 year warranty. These warranties exclude wheels and cables.



CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com





The MAXI+

Accessories



TRAY



OXYGEN TANK HOLDER



CANE HOLDER



SLOW DOWN BRAKE



ONE HAND BRAKE

Dolomite has a wide range of accessories to enhance the use of your walker. The Tray sits on the seat to give you a platform to carry food and beverages. The Oxygen Tank Holder replaces the basket and will safely support a B, C or D size tank. Attach your cane or crutch with the Cane Holder. Slow down brakes have a tension control knob to allow individual settings for each wheel. Slow down brakes are independent of the DABS braking system. One Hand brakes can control standard brakes from one handle. Replace the basket with a Sport-style Bag, it secures to the frame with Velcro and the large center compartment zips shut.

SPECIFICATIONS	LOW	STANDARD
Dist. Between Rear Wheels	22"	22"
Minimum Handle Height	28.5"	32.5"
Maximum Handle Height	35"	39.5"
Seat Height from Floor	22"	25.5"
Seat Dimensions	8.5" Width x 17" Height	
Overall Width	27.5"	27.5"
Overall Length	26"	29"
Width Between Handles	20.5"	20.5"
Folded Dimensions	33" Length x 27.5" Width x 11" Height	
Walker Weight (lbs.)	20.5	21
Weight Capacity (lbs.)	440	440
Product Number	D12122	D12120

CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com

