



THE ROHO® HYBRIDSELECT™ LOW AIR LOSS MATTRESS OVERLAY SYSTEM

SKIN PROTECTION

Introducing HybridSelect™ with HealFloat™! A new low air loss mattress overlay system that revolutionizes the overlay market by providing customized comfort and optimal skin protection. HybridSelect with HealFloat combines the power of SelectAir® low air loss systems with the ultimate skin protection of a ROHO® HIGH PROFILE® Cushion in the heel area to create an optimal environment for skin healing and protection. Simply remove the cushion to use on any chair for added skin protection.

VERSATILITY

HybridSelect is versatile. It sets up in minutes on top of any existing mattress – no need to store mattresses.

SIMPLICITY

This powerful low air loss system gives you advanced, easy-to-operate settings in a simple overlay. With the help of HybridSelect, it's never been this easy to protect patients' skin and keep them comfortable.



HYBRIDSELECT™ LOW AIR LOSS MATTRESS OVERLAY SYSTEM

SKIN PROTECTION.

The HybridSelect low air loss mattress overlay system creates an optimal environment for skin healing and protection. The zoned overlay utilizes an internal manifold which allows three independent mattress sections to be set using one button on the control unit. HybridSelect manages moisture and temperature in order to promote skin healing with "true" low air loss - allowing air to escape through the top of the mattress.

VERSATILITY.

Don't sacrifice care due to cost. Get the ultimate versatility in pressure distribution with the HealFloat feature. HealFloat provides optimal heel protection while the patient is lying in the bed. Simply remove the ROHO HIGH PROFILE cushion in the heel area and place it on a chair for continued, maximum skin protection.

SIMPLICITY.

The HybridSelect overlay sets up in minutes on top of an existing mattress eliminating the inconvenience of storing the mattress. This easy-to-use low air loss system delivers maximum performance with many advanced features and a powerful control unit.



APPLICATIONS

THE HYBRIDSELECT SYSTEM IS MOST EFFECTIVE FOR USERS WHO:

- require optimal pressure distribution.
- require increased heel protection.
- desire continued tissue protection while sitting in a chair.
- have minimal space; current mattress cannot be stored.

FEATURES & BENEFITS

POWERFUL:

Control unit provides constant inflation and can be adjusted to suit your needs.

ADJUSTABLE:

Mattress overlay and DRY FLOATATION® cushion can each be adjusted to fit any body shape.

COMFORTABLE:

Adjustable control unit settings for optimum comfort.

DURABLE:

High grade materials provide long product life even if patient's needs change.

WASHABLE:

Mattress overlay and cushion are easily cleaned. (See operations manual for specific cleaning instructions.)

SPECIFICATIONS

INCLUDES:

Mattress overlay, control unit, sheet, ROHO HIGH PROFILE cushion with cover, cushion pump, operations manual.

CONSTRUCTION:

Mattress overlay: Zoned overlay
 Overlay Sheet: Vapor permeable waterproof
 Control Unit: Variable speed blower motor
 Voltage: 120 AC, 50/60 Hz;
 Amps: Less than 5
 HIGH PROFILE Cushion: Neoprene

SIZES:

Mattress Overlay: 4" Height x 80" Length x 36" Width (10cm Height x 213cm Length x 91cm Width).
 Cushion: 18" x 16" (48cm x 41cm).

PERFORMANCE:

Pulsation: 1-minute, 5-minute intervals
 Mattress Overlay Weight Limit: 300 lbs / 136 kg.

WARRANTY:

1-Year Limited Warranty on Mattress System.
 2-Year Limited Warranty on ROHO Cushion.

MEDICARE HCPCS CODE:

E0372

We highly recommend distribution through a certified rehab professional when ordering and fitting a ROHO Group product.

MORE PRODUCTS FROM ROHO



CONTOUR SELECT™ CUSHION



HIGH PROFILE® QUADRO® SELECT® CUSHION



PRODIGY MATTRESS OVERLAY®

FOR ADDITIONAL INFORMATION, PLEASE CONTACT:



100 North Florida Avenue • Belleville, IL 62221-5429 • USA • www.therohogroup.com
 U.S. & Canada: 1 800 851 3449 • 1 618 277 9173 • fax 1 888 551 3449
 Outside the U.S. & Canada: 1 618 277 9150 • fax 1 618 277 6518

©2006 The ROHO Group, Inc. The following are registered trademarks of The ROHO Group: ROHO®, shape fitting technology®, HybridSelect™, HealFloat™, HIGH PROFILE®, SelectAir®, DRY FLOATATION®, QUADRO® SELECT®, and CONTOUR SELECT™.

The ROHO Group, Inc. has a policy of continuous product improvement and reserves the right to amend specifications presented in this sell sheet.