

WEIGHS ONLY 14.5 lbs.!



Fly-Lite Aluminum Transport Chair

Available
1 Pack

- DFL19-BL 19" Blue Frame and Blue Plaid Upholstery, (Has replaced TLC19-BL), 1/cs
- DFL19-BLK 19" Black Frame and Black/White/Red Plaid Upholstery, (Has replaced TLC19-BLK), 1/cs
- DFL19-RD 19" Burgundy Frame and Black Upholstery, (Has replaced TLC19-RD), 1/cs
- LELR-TF Swing-away Elevated Legrest, Tool Free, (Optional front riggings), 1pr/cs

FEATURES:

- Weighs only 14.5 lbs.
- Aluminum frame is lightweight and strong.
- Comes with carry pocket on back rest.
- Comes with seat belt for added safety.
- Attractive plaid style nylon upholstery.
- Tool free, height adjustable swing-away footrests come standard.
- 8" casters in rear. 6" casters in front.
- Back folds down for storage and transport and features deluxe back release.
- Comes with aluminum rear wheel locks.
- Weight Capacity: 300 lbs.
- Limited Lifetime Warranty.

CARRY BAG AVAILABLE
SEE PAGE 79

SPECIFICATIONS

| Item # | Overall Width (open) | Overall Width (closed) | Overall Height | Seat To-Floor Height | Back Height | Width Seat | Depth Seat | Overall Length w/ Riggings | Product Weight | Carton Shipping Weight | Weight Limit |
|---------------------|----------------------|------------------------|----------------|----------------------|-------------|------------|------------|----------------------------|----------------|------------------------|--------------|
| DFL19 SERIES | 22" | 9" | 36" | 18" | 18" | 19" | 15.5" | 33" | 14.5 lbs. | 19.5 lbs. | 300 lbs |

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

To order call toll free: 877.224.0946

WHEELCHAIRS

drive