




**Deluxe Fly-Weight Aluminum Transport Chair**

Available **1 Pack**

- FW17BL 17" Blue Frame and Black Upholstery, 1/cs
- FW17GR 17" Green Frame and Black Upholstery, 1/cs
- FW17SL 17" Silver Frame and Black Upholstery, 1/cs
- FW17DB 17" Black Frame and Tan Plaid Upholstery, 1/cs
- FW17BG 17" Green Frame and Burgundy Upholstery, 1/cs
- FW19BL 19" Blue Frame and Black Upholstery, 1/cs
- FW19GR 19" Green Frame and Black Upholstery, 1/cs
- FW19SL 19" Silver Frame and Black Upholstery, 1/cs
- FW19DB 19" Black Frame and Tan Plaid Upholstery, 1/cs
- FW19BG 19" Green Frame and Burgundy Upholstery, 1/cs
- LELR-TF Swing-away Elevated Legrest, Tool Free, (Optional front riggings), 1pr/cs

**FEATURES:**

- Weighs Only 19 lbs. about 30% less than traditional transport chairs.
- Aluminum frame for lightweight strength and durability.
- Composite wheels are lightweight and maintenance free.
- Deluxe back release.
- Deluxe all aluminum rear wheel locks.
- Comes with carry pocket on back of rest.
- Both 17" and 19" come with tool free, swing-away footrests.
- Comes with seat belt for added safety.
- Soft plush upholstery.
- Back folds down for easy storage and transport and features deluxe back release.
- 8" casters in front and rear.
- Weight Capacity: 300 lbs.
- Limited Lifetime Warranty.



**UPHOLSTERY**



**CARRY BAG AVAILABLE**  
SEE PAGE 79

**SPECIFICATIONS**

Item #	Overall Width (open)	Overall Width (closed)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length w/Riggings	Product Weight	Carton Shipping Weight	Weight Limit
<b>FW17 Series</b>	20"	9"	36"	19"	18"	17"	16"	33"	19 lbs.	24 lbs.	300 lbs.
<b>FW19 Series</b>	22"	9"	37.25"	19"	18"	19"	16"	33"	19 lbs.	24 lbs.	300 lbs.

Due to product improvements and changes, all images in this catalog are subject to change.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.