

WEIGHS ONLY 26 lbs.!



TR17

**DURABLE
CHROME FINISH!**



**CONVENIENT CUP HOLDER
CAN EASILY ATTACH TO CHAIR!
SOLD SEPARATELY MODEL # STD51040S**

**COMES WITH TOOL FREE
ADJUSTABLE FOOTRESTS**

Deluxe Go-Kart Steel Transport Chair (Chrome)

**Available
1 Pack**

TR17 17" Chrome Frame, Black Upholstery, 1/cs
TR19 19" Chrome Frame, Black Upholstery, 1/cs
LELR-TF Swing-away Elevated Legrest, Tool Free, 1pr/cs,
(Optional front riggings)

FEATURES:

- Weighs 26 lbs.
- Carbon steel frame with double coated chrome provides an attractive, chip-proof, easy to maintain finish.
- Composite 8" wheels are lightweight and maintenance free.
- Both 17" and 19" come with swing-away footrests.
- Padded armrests provide additional comfort.
- Nylon upholstery is comfortable, lightweight and easy to clean.
- Back folds down for storage and transport.
- Aluminum casted back release hinge allows back to fold down for easy storage and transport.
- Comes with seat belt for added safety.
- Comes with carry pocket on back rest.
- Comes with upholstery handles to fold chair.
- Standard rear wheel locks.
- Limited Lifetime Warranty.

SPECIFICATIONS

Item #	Overall Width (open)	Overall Width (closed)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length w/Rigging	Product Weight	Carton Shipping Weight	Weight Limit
TR17	20"	9"	36"	19"	18"	17"	16"	33"	26 lbs.	31.5 lbs.	300 lbs.
TR19	22"	9"	37.25"	19"	18"	19"	16"	33"	26 lbs.	31.5 lbs.	300 lbs.

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.