



## Super Light, Folding Transport Chair with Carry Bag and Flip Back Arms Available **1** Pack

SL18 Silver Frame and Blue Upholstery, 1/cs  
 835N Replacement Carry Bag for SL18  
 (Has replaced 839), 1/ea

### FEATURES:

- Padded flip back armrests.
- Weighs only 19 lbs. Easy to push or transport.
- Aluminum frame is lightweight and strong.
- Composite wheels are lightweight and maintenance free.
- Folds like a conventional folding chair for added convenience.
- Comes with carry bag (item #835N).
- Fold down footrests.
- Standard rear wheel locks.
- Attractive frame with durable, lightweight nylon upholstery.
- Weight Capacity: 250 lbs.
- Limited Lifetime Warranty.



TOOL FREE  
 REMOVABLE,  
 FRONT CASTERS



835N

FOLDS LIKE A CONVENTIONAL FOLDING CHAIR!



COMES WITH CARRY BAG WITH SHOULDER STRAP!



### SPECIFICATIONS

Item #	Overall Width (open)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length	Product Weight	Carton Shipping Weight	Weight Limit
SL18	26"	37"	20"	17"	18"	17"	30"	19 lbs.	24 lbs.	250 lbs.

Due to product improvements and changes, all images in this catalog are subject to change.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.