

# Trotter Mobility Chair



The Trotter is a lightweight portable positioning chair that provides transport solutions for parents and children. A wide range of options provide custom positioning for your individual needs and its simple adjustments allow the Trotter to grow with your child.

<b>Model #:</b>		<b>TR 1200</b>	Trotter, 12"
Sold 1 per box		<b>TR 1400</b>	Trotter, 14"
		<b>TR 1600</b>	Trotter, 16"
		<b>TR 1800</b>	Trotter, 18"

**Weight Capacity:** 75/100/175/250 lbs

**Warranty:** Limited Lifetime

**Accessories:** See page 32-33

## Features

- Lightweight folding frame with closure strap that holds the chair shut when folded.
- Black removable and washable fabric.
- Crash tested for transit. The Trotter has passed a dynamic crash test done with a 130 lb. front-facing occupant tested at 30 mph.
- Adjustable seat angle - 15° and 22°.
- Adjustable back angle - 85°, 90° and 95°.
- 8" front casters and 10" rear wheels with toggle lock.
- Adjustable seat depth.
- Curb assist lever facilitates curbside navigation.
- Positioning pelvic belt and five point harness with padded strap covers. The reversible padded strap covers increase comfort and one side of cover is imitation lamb's wool and the other is velour.
- Swing-away footrest are independently adjustable and removable.
- Limited Lifetime Warranty.



**CURB ASSIST LEVER**

## Specifications

Item #	Seat Width	Seat Depth	Back Height	Seat Seat-to Footplate	Seat to-Floor	Product Weight	Carton Shipping	Weight Limit Weight	Weight Limit w/Transit Option
TR 1200	11.75"	11.5" - 13.5"	20"	8" - 22"	22"	19 lbs.	28 lbs.	75 lbs.	75 lbs.
TR 1400	13.75"	13" - 15"	23"	9" - 22"	22"	21 lbs.	30 lbs.	100 lbs.	100 lbs.
TR 1600	15.5"	13.5" - 15.5"	24"	10" - 22.5"	22"	24.5 lbs.	33.5 lbs.	170 lbs.	175 lbs.
TR 1800	17.5"	15" - 17"	25"	11" - 23"	22.5"	30 lbs.	39 lbs.	250 lbs.	250 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

 = 1-Pack Available