





HCPCS
E1238

Wallaby Pediatric Folding Wheelchair

Item #	Product Seat Dimensions	Product Overall Width Open	Product Cane Height	Color	UOM
WB 1200	12" (W) x 12" (D)	20"	14.5"	Red 	1/ea
WB 1400	14" (W) x 12" (D)	22"	15"	Blue 	1/ea

specifications

Overall Width (Closed)	12"
Seat to Floor (Front)	19"
Seat to Floor (Rear)	18"
Seat to Footrest	10" - 15"
Push Handle Height	33" - 39"
Armrest Height	7"
Overall Length w/ Footrests	37"
Overall Length without Footrests	28.5"
Back Angle	90°
Product Weight w/ Footrests	37 lbs.
Product Weight w/out Footrests	33 lbs.
Weight Capacity	150 lbs.
Weight Capacity when Used in Bus Transport	80 lbs.
Warranty	Limited Lifetime

- WC19 Bus Transit Approved when used in conjunction with WB 8021 and WB 8030.
- Powder-coated colored steel frame with black cross brace.
- Easy-to-clean lightweight durable nylon upholstery.
- Height-adjustable push handles.
- Flip-back padded desk length arms.
- Pelvic belt with reversible cover.
- 22" rear and 6" front swivel solid rubber tires mounted on composite wheels.
- Push-to-lock brakes.
- Swing-away footrests with plastic plates and calf straps.