## FlexCHAIR<sup>TM</sup> Series







Pivoting arms with height adjustability



## FreedomFlex® Chair

- VBack Design® patent pending exclusive feature
  that keeps the resident centered in the chair while
  adding comfort and support while providing facilities
  the ability to standardize on a chair to accommodate
  more resident sizes
- Arm Rest height adjustable arm rest provide flexibility to customize height as needed and as an added benefit the pivoting motion allows an easier transfer
- Pressure Relief Strapping provide additional support and comfort by "flexing" to the individuals unique contour

- Steering Casters fixed rear casters allow the resident to self propel in a straight direction
- Adjust-ability customize the chair to the right height and seat depth based on the residents needs
- Independent Leg Rest optional bolt on accessory allows you to use any standard wheelchair leg rest (no costly proprietary design)
- Made in the USA an affordable quality product with a commitment to keeping inventory and parts on hand.(no customs or long lead times)



## Freedom Flex CHAIR SPECIFICATIONS

Frame .......Gauge steel 11ga - 16ga

Degree of Tilt......18°

Back Seat Height ......32"

Weight Capacity......350 LBS.

Warranty......3 year mechanical frame 1 year on components



## FlexChair Options



**FlexTilt**®

PROUDLY MADE IN THE U.S.A.



**FreedomFlex**®



Flex Rocker®

Coming soon in 16", 20", 24"



FALLS PREVENTION



CARE

is at the forefront of innovating medical technology for post acute care. We are focused on partnering with you to help address healthcare challenges in the following key areas—

MED MIZER®







REHABILITATION

WOUND CARE