

**VICTORY<sup>®</sup>**  
**ES 9**

ECONOMICALLY SOUND SCOOTER  
WITH IMPRESSIVE FEATURES  
AND PERFORMANCE

## VICTORY<sup>®</sup> ES 9

The **Victory<sup>®</sup> ES 9** is designed for durability at a maximum value. Features include an adjustable tiller, black non-scuffing tires, an LED headlight, color-through shrouds, rugged front plastics, and a crushproof basket.



**Pride<sup>®</sup>**  
Mobility Products Corp.  
— LIVE YOUR BEST —

800-800-8586 (US) • 888-570-1113 (Canada) • [www.pridemobility.com](http://www.pridemobility.com)

# Victory® ES 9 3-wheel

Model number:	S93
Weight capacity:	300 lbs.
Maximum speed <sup>1</sup> :	Up to 5 mph
Ground clearance <sup>2</sup> :	2.5"
Overall length <sup>2</sup> :	43"
Overall width <sup>2</sup> :	22.5"
Turning radius:	39.25"
Range per charge <sup>1,5</sup> :	11 miles
Tire size:	Front: 8" solid Rear: 9.5" solid
Total weight without batteries <sup>2,6</sup> :	101 lbs.
Heaviest piece when disassembled <sup>2</sup> :	Front: 45.5 lbs.
Standard seating:	Type: Foldable molded plastic Weight: 21.5 lbs. Material: Black vinyl Dimensions: 18" overall width 17" effective depth
Seat-to-deck range:	17.25" - 18.75"
Drivetrain:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual braking system:	Regenerative and electromechanical
Battery requirements <sup>3,4</sup> :	Two 12-volt, deep cycle Size: 18 AH Weight: 15 lbs. each
Battery charger:	3-amp, off-board

Colors:



Red



Blue

<sup>1</sup>Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.

<sup>2</sup>Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%.

<sup>3</sup>AGM or gel cell type required.

<sup>4</sup>Battery weight may vary based on manufacturer.

<sup>5</sup>Tested in accordance with ANSI/RENSA, WC Vol 2, Section 4, and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

<sup>6</sup>Includes standard seat weight.

## Features

- Pride's exclusive stylish, lightweight, low-profile, non-scuffing wheels
- Disassembles for transport or storage
- Bright LED headlight
- Seat features pinchless hinges and eight-position seat rotation
- Adjustable tiller
- 300 lbs. weight capacity

## Options


- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Weather cover
- Saddlebag

## Pride™ Lifts & Ramps

- Backpacker® Plus
- Outlander
- SilverBoom 250
- Backpacker® MV



Outlander



Backpacker® Plus

These lifts are compatible with the Victory® ES 9.