




3 Panel Privacy Screen

13508 1/cs

FEATURES:

- Sturdy, 1" diameter anodized aluminum tubing.
- Screens are 6 mil. flame retardant bacteriostatic white vinyl.
- Four, 3" hooded casters.
- Easy to maintain and clean.
- Overall Height: 70 1/2".
- Product Weight: 15 lbs.
- Shipping Weight: 20 lbs.
- Each hinged panel measures: 56" (H), 27 1/2" (W)
- Limited Lifetime Warranty.



Two-Piece Foam Cervical Collar (Philadelphia® Style)

3005-SM Small, 1/ea
 3005-MD Medium, 1/ea
 3005-LG Large, 1/ea
 3005-XL Extra Large, 1/ea

FEATURES:

- High density foam with Velcro®-type closure provides maximum support.
- Restricts movement of head to immobilize the cervical spine.
- Oversized trachea opening for emergency tracheotomies or carotid pulse monitoring.
- Molded design provides maximum support and patient comfort.
- Quick, easy application.
- Water-resistant.
- Easy to clean with mild detergent.
- Two-piece design allows the collar to be stored flat.
- 90 Day Limited Warranty.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.