FlexiFit*407

PREMIUM NASAL MASK

Specification Sheet

A Nasal Delivery Mask






Incorporating FlexiFit™
Technology with auto-contouring
properties, the FlexiFit™ 407
Nasal Mask comes with a unique
Glider™ Strap mechanism for
maximum comfort and stability.
Other unique features ensure a
premium sleep experience.





Nasal Delivery - *FlexiFit*™407




The FlexiFitTM407 Premium Nasal Mask offers an over-the-nose mask solution that fits better, feels lighter and is easier to use. The StretchgearTM Headgear provides even greater comfort and a maintenance-free bias diffuser ensures quieter operation. The FlexiFitTM407 Nasal Mask comes as an Out-of-Box Solution fully assembled and ready for use.

UNIQUE MASK FEATURES FOR GREATER FREEDOM, COMFORT AND COMPLIANCE



FlexiFit™ Technology Auto-contouring fit

With Foam Cushion
FlexiFit™ Technology
automatically adjusts to a wide
range of nasal and facial contours
with simple one-step ease
of fitting and no complicated
adjustments.



Glider™ Strap

The sliding Glider™ Strap provides freedom of movement while maintaining a seal. Mask force remains equal when the head rotates from side to side, minimizing the occurrence of leaks.



Stetchgear™ Headgear Deluxe mask suspension system

This is a special elastic headgear segment that maintains an equal mask force when the head tilts up or down.



Premium Frosted Silicone Seal Optimizing seal technology

This seal technology has a soft, smooth feel to slide onto the face, and its frosted texture optimizes the mask seal.



Efficient Bias Diffuser Quiet usage

Maintenance Free
This advanced air diffuser
system prevents air from
blowing onto the sleeping
partner and decreases noise that
could potentially disrupt sleep.



Glider™ Clip One-step remova

Simple one-step removal of headgear. No ongoing headgear adjustments required.



Hidden Strap Release One-step removal

Easy one-step removal of headgear. Enables reassembly without having to make fitting adjustments.

CONTENTS

- Complete mask with headgear assembled and ready for use
- Single one size Silicone Seal designed to fit most patients
- User friendly Fitting Guide

REORDER CODE: HC407

FOR MORE INFORMATION, PLEASE CONTACT YOUR LOCAL FISHER AND PAYKEL HEALTHCARE REPRESENTATIVE:

 $\textbf{International} \hspace{0.2cm} \textbf{PO Box 14 348}, \textbf{Panmure}, \textbf{Auckland 1134}, \textbf{New Zealand Tel: +64 9 574 0100} \hspace{0.2cm} \textbf{Fax: +64 9 574 0158} \hspace{0.2cm} \textbf{Email: info@fphcare.com Web: www.fphcare.com Web: www.fphcare.com$

USA Tel: 1800 446 3908 +1 949 470 3900 Fax: +1 949 470 3933 **Australia** Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 78 644 0924 Fax: +31 78 644 0914 **China** Tel: +86 20 8363 4286 Fax: +86 20 8363 3595 **France** Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 **Germany** Tel: +49 7182 93777 0 Fax: +49 7182 93777 99 **India** Tel: +91 80 5123 6041 Fax: +91 80 5123 6044 **Irish Republic** Tel: 1800 409 011 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 3661 7205 Fax: +81 3 3661 7206 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **UK (EU Authorised Representative)** Tel: +44 1628 626 136 Fax: +44 1628 626 146 **Nthn Ireland** Tel: 0800 132 189