FlexiFit¹432

FULL FACE MASK



Specification Sheet

A Full Face Delivery Mask






Our second generation Full Face Mask incorporates a soft Foam Cushion that works together with the highly contoured Silicone Seal to provide improved sealing technology and a more comfortable fit.





Full Face Delivery - *FlexiFit*™432




The FlexFitTM 432 Full Face Mask offers an improved seal and more comfortable fit with the inclusion of a soft Foam Cushion that interacts with the Silicone Seal. This alternative to nasal or oral delivery comes fully assembled and ready for use.

UNIQUE MASK FEATURES FOR GREATER FREEDOM, COMFORT AND COMPLIANCE



FlexiFit™ Technology

With Foam Cushion FlexiFit™ Technology automatically adjusts to a wide range of nasal and facial contours with simple one-step ease of fitting and no complicated adjustments. The Foam Cushion improves the sealing technology and comfort.



Glider™ Strap Anti-leak system

The sliding Glider™ Strap provides freedom of movement while maintaining a seal.

Mask force remains equal when the head rotates from side to side, minimizing the occurrence of leaks.



Jnder-Chin Design Extra stability

An Under-Chin fit provides stability by providing a constant point of reference and maintains the seal even with movement of the jaw during the night.



Efficient Bias Diffuser

Maintenance Free
This advanced air diffuser
system prevents air from
blowing onto the sleeping
partner and decreases noise that
could potentially disrupt sleep.



Premium Frosted Silicone Seal Optimizing seal technology

This seal technology has a soft, smooth feel to slide onto the face, and its frosted texture optimizes the mask seal.



Crown Strap Extra stability

Complements the Under-Chin Design to maximize stability and eliminate the need for a chin strap.



Easy-Release Strap One-step removal

Simply pull the Easy-Release Cord for one-step removal of headgear. No ongoing headgear adjustments required.

CONTENTS

- Complete mask with headgear assembled and ready for use
- The FlexiFit[™] 432 Full Face comes in three separate sizes, Small, Medium and Large
- User friendly Fitting Guide

REORDER CODE: HC432

FOR MORE INFORMATION, PLEASE CONTACT YOUR LOCAL FISHER AND PAYKEL HEALTHCARE REPRESENTATIVE: