

Freedom to move

Opus.360

NASAL PILLOWS MASK

The Opus^{*360} Nasal Pillows Mask is small, quiet and superbly comfortable, offering unprecedented freedom. The slim-line design allows for almost any sleeping position without mask displacement. Its innovative 360-degree pivoting ball-and-socket elbow provides greater versatility. With its new Silicone Pillows and facially Contoured Mask Frame for extra stability, patient comfort is paramount.

Elegant in design and simplicity, let the Opus⁻³⁶⁰ Nasal Pillows Mask from Fisher & Paykel Healthcare work for you today!


NEW:: Contoured Mask Frame

Facially contoured, the three-dimensional support structure conforms naturally to provide maximum stability and comfort.



Fisher & Paykel HEALTHCARE

Direct Nasal Delivery – Opus



The Opus^{**360} Nasal Pillows Mask is small, quiet and superbly comfortable, with anatomical pillow-style Silicone Seals and a facially Contoured Mask Frame. A ball-and-socket elbow combines with an innovative anchoring system to give unprecedented freedom of movement and stability. The Opus^{**360} Nasal Pillows Mask is a complete, Out-of-Box Solution, fully assembled and ready for use with alternative sizes of Silicone Pillows included.

UNIQUE MASK FEATURES FOR GREATER FREEDOM, COMFORT AND COMPLIANCE



NEW:: Anatomically Shaped Pillows Mask comfort

The Anatomically Shaped Pillows provide easy fitting with an excellent seal. The contours of the nasal pillows ensure greater comfort. The range of sizes provided (S,M,L) ensures comfort and seal are not compromised.



NEW:: Contoured Mask Frame Mask stability

The facially Contoured Mask Frame, with three-dimensional support structures, conforms naturally to the shape of the face to provide a comfortable fit and reliable support. Its slim-line design allows side-sleeping without mask displacement.



Tube Anchoring System Mask stability

The innovative pivot system allows the tube to be anchored at the top or side, or left free to track. Elegant in design and simplicity, it firmly holds the tube in the desired position. Alternatively, the tube can be left unanchored so that the balland-socket elbow allows true 360-degree freedom.



Efficient Exhaus

The advanced air diffuser system is quiet, minimizing the noise and draft of the exhaust flow to prevent both sleeper and partner disturbance.



Soft-Lip™ Cushion Mask comfort

The Silicone Seal is soft on the nose and supple on the upper lip. The built-in Soft-Lip™ Cushion conforms naturally to the shape of the upper lip, for maximum comfort.



Premium Frosted Silicone Seal Optimizing seal technology

The enhanced seal technology has a soft, smooth feel and slides easily into position. Its frosted texture optimizes mask seal performance.



Hidden Strap Release
One-step removal

Easy one-step removal of Headgear for cleaning. Enables reassembly without having to make readjustments. Simply slide the Side Strap down and out of the Hidden Strap Release hooks.

CONTENTS

- Complete Mask with Mask Frame assembled and ready for use
- All three sizes of Silicone Seal included (S,M,L)
- User-friendly Fitting Guide

REORDER CODE: HC482

FOR MORE INFORMATION, PLEASE CONTACT YOUR LOCAL FISHER & PAYKEL HEALTHCARE REPRESENTATIVE:

International P 0 Box 14 348, Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.com Web: www.fphcare.com Australia Tel: +61 3 9879 5022 Fax: +61 3 9879 5022 Fax: +61 3 9879 5222 Austria Tel: 0800 29 31 23 Fax: 0800 29 31 22 Benelux Tel: +31 78 644 0924 Fax: +31 78 644 0924 Fax: