


# Aqua Creek Products

The Total Solution for Pool and Spa Access



## Item# F-12PPL-HD-AT1

The Patriot™ Lift is the only completely portable aquatic lift that is independently verified to meet ADA guidelines. It has a lifting capacity of 450lb\* (204kg). Our fixing kit (included) provides a solution to the revisions made by the Department of Justice.

This lift can easily be removed from poolside when not needed. It features a compact, low profile frame that does not intrude into the pool when not being used. The unique sliding mechanism allows the lift to be deployed and returned to the cart without fastening hardware.

U.S. Patent No.  
7,249,386 B2

U.S. Patent No.  
5,790,995

U.S. Patent No.  
D507,769 S



Push lift to pool side, set brakes, slide into place, and it's ready to use.

### OPTIONS

- **Upgrade Pack - Item #F-12PPAP**  
Includes headrest, chest strap, battery & cover
- **Headrest - Item #F-422PLH**
- **Chest strap - Item #F-423CS**
- **Lockable cover - Item #F-12PPC**
- **Pull-out leg rest - Item #F-105LAR**
- **Cycle Attachment - Item #F-019CA**



### FEATURES

- 450 lb\* (204 kg) weight capacity
- Independently verified to ADA standards
- Completely portable unit- no anchor required
- Rechargeable battery with wall mount charger
- Adjustable footrest & lap belt
- Dual flip-up arm rests
- Submersible remote with push button operation
- Limited lifetime structural warranty
- 5 year pro-rated electronics warranty

\*450 lb (204 kg) weight capacity for water drafts up to 6"  
400 lb (181 kg) capacity for other applications.



**PROUDLY SOLD BY:**