

$Versamode^{TM} - Drop Arm$

- Provides transport, transfer, and commode functions all in one unit
- Arms drop independently below seat level for safe lateral transfers
- Cushioned armrests enhance user comfort
- Swing-away footrests can be extended or detached
- 5" rear swivel, locking casters facilitate maneuverability in tight spaces
- · Padded seat cover also serves as a backrest
- Maximum Weight Capacity 300 lb EVENLY DISTRIBUTED
- Comes unassembled
- Designed to meet the requirements of HCPCS code: E0165
- · Limited lifetime warranty

6810A 1 ea





Specifications