

# Aqua Creek Products

The Total Solution for Pool and Spa Access

Pool Access Chairs are an ideal solution for swimming pools with ramps or zero depth entries. These chairs can be used in showers, locker rooms and pool environments. The Pool Access Chairs feature an anti-tip design with non-rusting parts constructed from medical grade PVC. Swing out armrest and retractable footrests make lateral transfers easy for any patron. Chairs ship truck freight fully assembled.

## Features:

- **250 lb (113 kg) weight capacity - Pool Access Chair**
- **400 lb (181 kg) weight capacity - Heavy Duty Chair**
- **Durable non-corrosive PVC construction**
- **Stainless steel hardware**
- **Plastic rear wheels with rubber tread**
- **Anti-tip design**
- **Swing out armrest for lateral transfers**
- **One year warranty**

## Options:

- **Pool Access Chair with Optional Mesh Seat**  
*Item #F-011PAC-MS*
- **Chest strap -** *Item #F-426CS*
- **Seat belt -** *Item #F-11PACSB*



## Pool Access Chair

*Item #F-011PAC*



## Heavy Duty Pool Access Chair

*Item #F-400PWC*

**PROUDLY SOLD BY:**